

## OUR IMPLANT SYSTEMS

RELIABLE, EFFICIENT, SAFE





**iRES**<sup>+</sup>

Swiss  
implant  
quality







## INDEX

### Implant systems

ShapeOne®	p. 04
iMAX®	p. 05
iMAX® NHSIC	p. 06
iMAXMUA®	p. 07
Volution®	p. 08
ShapeMini	p. 09
Price list	p. 10
Colors implant systems	p. 11
Implant connections	p. 14
Management of soft tissues	p. 15
Surface treatment	p. 16
Cytotoxicity test	p. 17
Cold plasma decontamination	p. 17
Sterilization & packaging	p. 17

### Prosthetic solutions

Internal hex	p. 20
External hex	p. 28
Conical connection	p. 34
Internal octagon	p. 38
ShapeMini	p. 40
MUA components	p. 41
Accessories, rings and retaining caps	p. 43
Digital prosthetic solutions	p. 45

### Kit and instruments

Surgical kit	p. 49
Drills	p. 50
Countersinks	p. 51
Taps	p. 51
Stops	p. 51
Dynamometer ratchet	p. 52
Surgical kit composition	p. 53
Prosthetic kit composition	p. 59
Ergo kit composition	p. 60
Sinus lift kit composition	p. 61

### Surgical protocol

Surgical protocol	p. 62
Surgical protocol for short implants	p. 69
Minimum implants size	p. 75


### Package leaflets

Implant systems	p. 77
Prosthetic components	p. 78
Rotary instruments	p. 79

### Certification

	p. 80
--	-------









## OUR IMPLANT SYSTEMS

---

Reliable  
Efficient  
Safe

SHAPEONE<sup>®</sup>

iMAX<sup>®</sup>

iMAX<sup>®</sup> NASIC

iMAXMUA<sup>®</sup>


VOLUTION<sup>®</sup>

SHAPEMINI<sup>®</sup>



Code	Diameter	Thread	Connection	Platform	Hex	Treatment
B	3.75 - 4.1 - 4.7	Triple	Internal hex	3.5 mm	2.5 mm	C - HYHA
Tn	3.75 - 4.1 - 4.7	Triple	Internal hex	3.5 mm	2.5 mm	C
T	3.75 - 4.1 - 4.7	Triple	Internal octagon	4.8 mm	3.1 mm	C

## EXCELLENT PRIMARY STABILITY ALSO IN BONE d4



1. Platform switching.
2. 1 mm machined neck to protect from bacterial attacks.
3. Tapered self tapping with cutting coronal thread for a better force discharge on cortical bone to provide excellent stability even in few millimeters.
4. 55° triple thread over the entire body implant with a pitch of 1,8 mm 4 (0,6 mm/thread); each turn allows to go down of 1,8mm, speeding up the insertion phase.
5. 3 apical aggressive cuts provide a better primary stability and centering of the implant and the possibility to change direction during its insertion.
6. Apex aggressive but rounded to protect the Schneider's membrane.

## TYPE OF SURFACE TREATMENT and CONNECTIONS



S1B - C  
Bone level



S1B - HYHA  
Bone level



S1TN  
Bone level o Tissue level




S1T  
Tissue level

Ø	ShapeOne B	ShapeOne T	ShapeOne Tn
Ø 3.75	h 8   10   11.5   13   16	h 8   10   11.5   13   16	h 8   10   11.5   13   16
Ø 4.1	h 6.5   8   10   11.5   13   16	h 4.5   6.5   8   10   11.5   13   16	h 4.5   6.5   8   10   11.5   13   16
Ø 4.7	h 6.5   8   10   11.5   13   16	h 4.5   6.5   8   10   11.5   13   16	h 4.5   6.5   8   10   11.5   13   16



Code	Diameter	Thread	Connection	Platform	Hex	Treatment
NHSI	3.3	Fine double thread	Internal hex	3.2 mm	2.1 mm	C - HYHA
NHSI	3.75 - 4.1 - 4.7 - 5.2	Fine double thread	Internal hex	3.5 mm	2.5 mm	C - HYHA
NHSE	3.3	Fine double thread	External hex	3.5 mm	2.4 mm	C - HYHA
NHSE	3.75 - 4.1 - 4.7 - 5.2	Fine double thread	External hex	4.1 mm	2.7 mm	C - HYHA

UNIVERSAL SYSTEM FOR ALL TYPES OF BONES



- 1. Platform switching.
- 2. 1 mm machined neck.
- 3. Cylindrical body.
- 4. 55° double thread over the entire body implant with a pitch of 1,2 mm (0,6 mm/thread); each turn allows to go down of 1,2mm, speeding up the insertion phase.
- 5. 2 apical cuts helicoidal.
- 6. Conical apex flat tip.

TYPE OF SURFACE TREATMENT and CONNECTIONS



NHSI - C  
Bone level



NHSI - HYHA  
Bone level



NHSE - C  
Bone level



NHSE - HYHA  
Bone level

Ø	NHSI	NHSE
Ø 3.3	h 10   11.5   13   16	h 10   11.5   13   16
Ø 3.75	h 8   10   11.5   13   16	h 6.5   8   10   11.5   13   16
Ø 4.1	h 6.5   8   10   11.5   13   16	h 6.5   8   10   11.5   13   16
Ø 4.7	h 6.5   8   10   11.5   13   16	h 6.5   8   10   11.5   13   16
Ø 5.2	h 6.5   8   10   11.5   13	h 6.5   8   10   11.5   13

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps

Digital prosthetic  
solutions

Drills, countersinks  
taps, stops

Surgical  
kit


Surgical  
protocol

Package  
leaflets



Code	Diameter	Thread	Connection	Platform	Hex	Treatment
NHSIC	3.3 - 3.75	Fine double thread	Cono-morse	Narrow	2.1 mm	C
NHSIC	4.1 - 4.7 - 5.2	Fine double thread	Cono-morse	Regular	2.5 mm	C

THE MORE RELIABLE CONICAL CONNECTION WITH DIFFERENT EMERGENCE PROFILES



- Platform switching, unique prosthesis for all the diameters.
- 1 mm machined neck.
- 55° double thread over the entire body implant with 5 a pitch of 1,2 mm (0,6 mm/thread); each turn allows to go down of 1,2mm, speeding up the insertion phase.
- 2 apical cuts helicoidal.
- Conical apex flat tips.
- 5° cone inclination on implant and abutment.

Range of heights, from 2.5 to 4 mm, healing srews, temporary transfers and permanent abutments

It's proven that the cone-morse connection creates smaller slits (1µm) of bacteria (1,1-1,5µm length, 2-6µm diameter). Cone-morse connection absorbs vibration and chewing stress by eliminating the unscrewing of the screws (0.37%).

Healing abutment with passing screw doesn't ruin the connection, like happens in traditional spin. Screws with cone-morse connection, block the microleakage during the healing period. Also straight MUA and ball abutments are provided with the passing screw.

TYPE OF SURFACE TREATMENT and CONNECTIONS




Ø	NHSIC
Ø 3.3	h 10   11.5   13   16
Ø 3.75	h 8   10   11.5   13   16
Ø 4.1	h 6.5   8   10   11.5   13   16
Ø 4.7	h 6.5   8   10   11.5   13
Ø 5.2	h 6.5   8   10   11.5   13



Code	Diameter	Thread	Connection	Platform	Angle	Treatment
NHSM 00	3.3 - 3.75 - 4.1	Double	One Piece	4.3 mm	0°	HYHA
NHSM 18	3.75 - 4.1	Double	One Piece	4.3 mm	18°	HYHA
NHSM 30	3.75 - 4.1	Double	One Piece	4.3 mm	30°	HYHA

ONE PIECE IMPLANT FOR IMMEDIATE LOADING



- 1. Hole for retaining screw.
- 2. 1.5mm neck.
- 3. Cylindrical body.
- 4. 55° double thread over the entire body implant with a pitch of 1,2 mm 3 (0,6 mm/thread); each turn allows to go down of 1,2mm, speeding up the insertion phase.
- 5. 2 apical helicoidal cuts.
- 6. Conical flat tip apex.

The use of the retention screw for MUA components, from the current 1.4 mm to 1.72 mm to eliminate the unscrewing possibility of MUA components.

Smaller size both of the cone and the supporting surface, from 4.8 mm to 4 mm. Reduced footprint for the prosthesis benefit.

TYPE OF SURFACE TREATMENT and CONNECTIONS



iMAX - NHSM00  
Tissue level



iMAX - NHSM18  
Tissue level




iMAX - NHSM30  
Tissue level

Ø	NHSM00	NHSM18	NHSM30
Ø 3.3	h 4.5   8   10   11.5   13		
Ø 3.75	h 4.5   8   10   11.5   13		h 11.5   13   16
Ø 4.1	h 4.5   8   10   11.5   13		h 11.5   13   16



Code	Diameter	Thread	Connection	Platform	Hex	Treatment
SVB	3.3	Double	Internal hex	3.2 mm	2.1 mm	C
SVB	3.75 - 4.1 - 4.7 - 5.2	Double	Internal hex	3.5 mm	2.5 mm	C

LARGE THREAD IMPLANT FOR A BETTER INSERTION FACILITY



1. Hexagonal internal connection with Friction Fit to avoid micromovements and bacterial infiltrations.
2. Machined implant neck with double platform switching and profile with toothed hermetic cover. It distributes the pressure on the cortical bone and ensures a seal that prevents the penetration of bacteria and the onset of peri-implantitis.
3. Implant with ogival design (bullet type). It ensures a higher penetration capacity.
4. Pronounced spires with cutting profile. They provide higher primary stability and self-tapping properties. Higher stability in soft bone 3/4 and less compression in hard bone 1/2).
5. Spiral counter-unload furrows. They increase primary 5 stability and prevent the progression of peri-implantitis, interrupting the spires that act as a chute for bacteria.
6. Tip with self-perforating auger design with grinding 6 effect. It produces bone shavings and stimulates osteogenesis.

TYPE OF SURFACE TREATMENT and CONNECTIONS




Volution - SVB  
Bone level

Ø	Volution SVB
Ø 3.3	h 8   10   11.5   13   16
Ø 3.75	h 6.5   8   10   11.5   13   16
Ø 4.1	h 6.5   8   10   11.5   13   16
Ø 4.7	h 6.5   8   10   11.5   13   16
Ø 5.2	h 6.5   8   10   11.5   13



Code	Diameter	Thread	Connection
SM	2.7	Mono	With ball connection in TIN
SMF	2.7	Mono	With fixed abutment

SELF-TAPPING IMPLANT WITH SINGLE THREAD



- 1. Ball attachment system (SM version), ideal for overdenture
- 2. Connection suitable for the cementation of the abutments (SMF version)
- 3. Self-tapping cylindrical body with single thread, pitch of 1.2 mm
- 4. Thread with lower angle of 43° and upper of 23°

TYPE OF SURFACE TREATMENT and CONNECTIONS



ShapeMini - **SM**  
With ball connection in TIN



ShapeMini - **SMF**  
With fixed abutment


Ø	SM	SMF
Ø 2.7	h 8   10   11.5   13   16	h 8   10   11.5   13   16




IMPLANT SYSTEMS						price list
	Implant system	Type	Ø	Connection	Treatment	Price
Implant systems	ShapeOne	B	3.75 - 4.1 - 4.7	Internal hex	C	130, <sup>00</sup>
Internal hex prosthetic solutions	ShapeOne	B	3.75 - 4.1 - 4.7	Internal hex	HYHA	165, <sup>00</sup>
External hex prosthetic solutions	ShapeOne	Tn	3.75 - 4.1 - 4.7	Internal hex	C	135, <sup>00</sup>
Conical connection prosthetic solutions	ShapeOne	T	3.75 - 4.1 - 4.7	Internal octagon	C	150, <sup>00</sup>
Internal octagon prosthetic solutions	iMAX	NHSI	3.3	Internal hex	C	135, <sup>00</sup>
	iMAX	NHSI	3.3	Internal hex	HYHA	165, <sup>00</sup>
	iMAX	NHSI	3.75 - 4.1 - 4.7 - 5.2	Internal hex	C	130, <sup>00</sup>
Shapemini prosthetic solutions	iMAX	NHSI	3.75 - 4.1 - 4.7 - 5.2	Internal hex	HYHA	165, <sup>00</sup>
MUA and OnePiece components	iMAX	NHSE	3.3	External hex	C	135, <sup>00</sup>
	iMAX	NHSE	3.3	External hex	HYHA	165, <sup>00</sup>
Accessories, rings and retaining caps	iMAX	NHSE	3.75 - 4.1 - 4.7 - 5.2	External hex	C	130, <sup>00</sup>
	iMAX	NHSE	3.75 - 4.1 - 4.7 - 5.2	External hex	HYHA	165, <sup>00</sup>
Digital prosthetic solutions	iMAX	NHSIC	3.3 - 3.75	Narrow	C	155, <sup>00</sup>
	iMAX	NHSIC	4.1 - 4.7 - 5.2	Regular	C	155, <sup>00</sup>
	iMAX MUA	MHSM 00	3.3 - 3.75 - 4.1	OnePiece	HYHA	195, <sup>00</sup>
Drills, countersinks taps, stops	iMAX MUA	MHSM 18	3.75 - 4.1	OnePiece	HYHA	225, <sup>00</sup>
	iMAX MUA	MHSM 30	3.75 - 4.1	OnePiece	HYHA	225, <sup>00</sup>
Surgical kit						
Surgical protocol	Volution	SVB	3.3 - 3.75 - 4.1 - 4.7 - 5.2	Internal hex	C	135, <sup>00</sup>
	ShapeMini	SM	2.7	Ball connection		55, <sup>00</sup>
Package leaflets	ShapeMini	SMF	2.7	With fixed abutment		55, <sup>00</sup>



# COLORS OF IMPLANT SYSTEMS AND SURFACE TREATMENTS




## LEGEND REF CODE







Swiss  
implant  
quality








## MORE RELIABLE CONNECTIONS



Prosthetic components with **Friction Fit** connection and one conical degree on the wall of the hex abutment, have been developed for **S1B**, **S1TN**, **SVB** and **NHSI** implant systems with internal hex connection. This connection ensures a “cold fusion” between implant and abutment if the retaining screw has been tightened at 30 Ncm. It eliminates micro movements and reduces bacterial infiltration

between implant and abutment. Prosthetic components with tilted cone to 16° fit perfectly with our **S1T** implant system. A “flat to flat” connection for external hex is especially suitable for **NHSE** implant systems iRES® offers different implant products to cover a wide range of treatment options and fulfill every surgical need.



The **NHSIC conical connection** has a 5° angle on the abutment and on the implant, and an emergence profile for the mucous attack. The cone-morse connection creates fissures (1µm) smaller than bacteria, absorbs vibration and chewing stress with the result that it eliminates the unscrewing of the

screws. There is no connection between implant and abutment in **NHSM ONEPIECE** system and this allows to completely eliminate bacterial infiltration. OnePiece connection is suitable for immediate loading.



MANAGEMENT OF SOFT TISSUES

The connection between implant and abutment is present on the entire implant line with the configuration to be finished as shown in figure 1 or classic as in figure 2, depending on the diameter of the flaring required, respectively 3.5 mm or 5 mm. An emergence profile is also available, figure 3, which ensures a more rounded connection

without sharp corners, creating a greater space for mucosa anchoring. This profile which guarantees conditioning of the soft tissues is contemplated on the whole prosthetic range, from healing screws right up to the final abutments. In figure 4 we can see the new econometric connection with 5° angle on the implant and on the abutment.



(figure 1)



(figure 2)



(figure 3)



(figure 4)

Implant systems
Internal hex prosthetic solutions
External hex prosthetic solutions
Conical connection prosthetic solutions
Internal octagon prosthetic solutions
Shapemini prosthetic solutions
MUA and OnePiece components
Accessories, rings and retaining caps
Digital prosthetic solutions
Drills, countersinks, taps, stops
Surgical kit
Surgical protocol
Package leaflets




	SURFACE TREATMENT		SLA type
Implant systems	<p>There is a relevant scientific literature* on how surface roughness characteristics influence cell behaviour. Compared to a smooth surface, topographical patterns smaller in size than a fibroblast cell (micro and nano topography) orient the arrangement of the cells and stimulate osteoblastic and platelet activity, accelerating the production of extracellular matrix and bone regeneration, and therefore the osseointegration of the dental implant. The three fundamentals of surface treatment of dental implants from a biological point of view are:</p> <ol style="list-style-type: none"> <li>1] control of surface topography to stimulate cellular response in an osteogenic direction;</li> <li>2] control of the chemical composition of the surface to promote cell colonization;</li> <li>3] control of biological contamination from adherent endotoxins so as not to interfere with the natural inflammatory response. For the surface treatment a sand-blasting process was used followed by a double acid attack. In the images, increasing the magnification, it can be seen how the macroscopic aspects of the screw (spire, cutting SLA surface treatment edge) are not affected by the treatment and that the surface is free from processing residue. The dual-beam roughness typical of SLA treatment can be clearly observed, which contains large cavities due to large grit blasting on which is superimposed the micro-roughness due to treatment with acids. The micro-roughness illustrated in the figures highlights the typical three-dimensional topography, which gives these surfaces “sponge-like” characteristics that are the basis of their excellent clinical performance. In fact, the very short peak-to-peak distance, about 1 micrometer, stimulates both the activity of osteogenic cells and the capillary penetration of the blood in the surface structure, offering very favorable characteristics to stimulate bone regeneration, as described in many articles on this topic. This unique combination of long-range roughness (large grit sand-blasting) and short-range (acid etching) is a substrate favorable to cell regrowth that adequately promotes cell differentiation. The level of roughness is <math>Ra\ 1.42 \pm 0.12</math>.</li> </ol>		
Internal hex prosthetic solutions			
External hex prosthetic solutions			
Conical connection prosthetic solutions			
Internal octagon prosthetic solutions			
Shapemini prosthetic solutions	<div>NECK MACHINED</div> <div>C</div> <p>Surface treatment on the body implant</p>		
MU and OnePiece components	<div>HYALURONIC ACID</div> <div>HYHA</div> <p>Partial surface treatment on the body implant with hyaluronic acid. Cold plasma decontamination*</p> <div>  <p>MAG 52 X WD 11.5 mm EHT 20.00 kV Signal A CZ BSD</p> </div> <div>  <p>MAG 200 X WD 11.0 mm EHT 20.00 kV Signal A CZ BSD</p> </div> <div>  <p>MAG 1.50 K X WD 11.5 mm EHT 20.00 kV Signal A SE1</p> </div> <div>  <p>Sa 0,50 µm overall mean value on a measuring area of 30x30 µm cold plasma decontamination</p> </div> <div>  <p>Sa 1,90 µm overall mean value on a measuring area of 30x30 µm sand-blasting, double etching, cold plasma decontamination</p> </div>		
Accessories, rings and retaining caps			
Digital prosthetic solutions			
Drills, countersinks taps, stops			
Surgical kit			
Surgical protocol			
Package leaflets	<p>* Valutazione della composizione chimica superficiale, della morfologia, della citotossicità e dell'adesione cellulare su impianti dentali. G. Cascardo, C. Cassinelli. Doctor OS 2005 Nov-Dic; 16 (9): 1091. Valutazione comparativa del trattamento di superficie in 5 sistemi implantari. M. Biasotto, M. Cadenaro et al. Università degli studi di Trieste. Quintessence International, Anno 18 - Maggio/Giugno 2002. RAPPORTO ISTISAN 01/15 - Valutazione del trattamento superficiale sulle prestazioni meccaniche a fatica di impianti in titanio plasma-sprayed e titanio sabbato e mordenzato. Rossella Bedini, Gior-gio de Angelis, Marco Tallarico, Rosario Ialpi, Umberto Romeo, Giuseppe di Cintio 2001, 33 p. RAPPORTO ISTISAN 08/32 - Valutazione microtomografica dell'area di possibile contatto osseo di sei tipologie diverse di impianti dentali. Rossella Bedini, Raffaella Pec-ci, Fabio Di Carlo, Alessandro Quaranta, Francesca Rizzo, Manlio Quaranta, G. Heimke, W. Schulte, B. d'Hoedt, P. Griss, C.M. Busing, D. Stock. The influence of fine surface structures on the osseointegration of implants. The International Journal of Artificial Organs 1982; 5(3): 207-212. Guy, M.J. McQuade, M.J. Scheidt, J.C. McPherson III, J.A. Rossmann, T.E. Van Dyke. In vitro attachment of human gingival fibroblasts to endosseous implant materials. Journal of Periodontology 1993 Jun; 64(6): 542-546.</p>		

CYTOTOXICITY TEST


EN ISO 10993-5: 2009, Valutazione Biologica dei Dispositivi Medici - Prove per la citotossicità in vitro

After treatment and decontamination, the implants proved to be perfectly cytocompatible, that is devoid of cytotoxic effects against L929 fibroblasts. In all wells, the cells always showed density and morphology fully comparable with those of the negative control. The fibroblasts proliferate homogeneously in contact with the implants as the


Material does not release any cytotoxic element. Moreover, multinucleated giant cells were never shown in significantly higher number than the negative control, indicating the absence of effects of an inflammatory type.



negative control



positive control



examined control

COLD PLASMA DECONTAMINATION

After the surface treatment, the implants are cleaned to remove processing residues by washing them with solvents and then subjecting them to a process of surface decontamination with cold plasma (Argon). The partially ionized Argon atoms (inert gas) act as an additional atomic sand-blasting that promotes the removal of organic contaminants and activates the ionization of surface atoms of titanium, improving the wettability of the implant. The treatment conditions adopted on shape1 implants offer the best characteristics considered important, according to the state of current knowledge\*, in the

processes of implant healing, both in terms of surface morphology and in terms of chemical composition (surface cleaning). Plasma cleaning, packaging in a controlled environment, the absolute respect of “clean” procedures, quality control tests of during the manufacturing process, play a fundamental role in the control of adherent endotoxins (biological cleaning), the main agent of immunological response to implant surfaces.

\* Valutazione del rapporto tra costo e qualità della pulizia superficiale di alcuni sistemi implantari in commercio Marco Morra, Clara Cassinelli, Giovanna Cascardo, Daniele Bollati, Nobil Bio Ricerche srl Via Valcastellana 26, 14037, Portacomaro (AT)

M. Morra, C.Cassinelli, Evaluation of Surface Contamination of Titanium Dental Implants by Lu-Sem: Comparison with XPS Measurements Surface and Interface Analysts, Vol. 25, 983-984 (1997).

STERILIZATION & PACKAGING

To preserve its integrity, the dental implant is housed in a vertical position inside a titanium cylinder anchored, by means of the closing cap, to the respective vial made of borosilicate glass for pharmaceutical use, complying with the European Pharmacopoeia in force. This vial really ensures the neutrality of the primary packaging due to the absence of release of contaminants during the sterilization phase. It is inserted in a blister of transparent polyglass sealed with heat-sealing lacquer-based Tyvek and packed in a cardboard box that also contains the instructions for use and the labels for the patient records, on which are printed the data that allow product traceability (code and batch number). All the product packaging Materials have

been tested, approved and certified.

Shape1® implants are supplied sterile, in a pack that allows their stability to be guaranteed for 5 years. The sterilization process is performed with gamma rays respecting the standards in force by qualified suppliers who use automated, safe and reliable systems, with continuous microbiological monitoring of the process.



\*European Pharmacopoeia, current edition, 3.2.1 Glass containers for pharmaceutical use.

Implant systems

Internal hex prosthetic solutions

External hex prosthetic solutions

Conical connection prosthetic solutions

Internal octagon prosthetic solutions

Shapemini prosthetic solutions

MUA and OnePiece components

Accessories, rings and retaining caps

Digital prosthetic solutions

Drills, countersinks, taps, stops

Surgical kit

Surgical protocol

Package leaflets









## PROSTHETIC SOLUTIONS

---

Surgical screws  
Prosthetic screws  
Healing screws  
Analog  
Scan abutment  
Indirect impression transfer  
Direct impression transfer  
Temporary straight abutments  
Abutments for welded technique  
Definitive straight abutments  
Sticking bases  
Angled abutments  
Anatomic angled abutments  
Castable abutments  
Straight mua abutments  
Angled mua abutments  
Ball abutments  
iRetor

SHAPEONE

iMAX

iMAX *NASIC*

iMAXMUA

VOLUTION





SHAPEMINI











--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Internal hex		
	iMAX NHSI 3.3	iMAX NHSI ShapeOne B Volution SVB
description   €		

<b>ANALOG</b>	Material TI-6Al-4V
---------------	--------------------

Analog 21.00 €	 S1BNIA	 S1BIA
Short analog for guided surgery 21.00 €	 S1BNIAS	 S1BIAS

<b>INDIRECT IMPRESSION TRANSFER</b>	*ITC and S1BRS1 screw included   **S1BRS1 screw included   Material TI-6Al-4V
-------------------------------------	---

Ø 3.5 mm 65.00 €	 S1BN135ITC *	 S1B135ITC *
Ø 5 mm 65.00 €		 S1B150ITC *
Indirect transfer cap 20.00 €	 ITC	
Multifunction abutment 55.00 €	 Ø 3.5 mm S1BN1A35 **	 Ø 5 mm S1B1A50 **
Multifunction abutment 55.00 €	 Ø 3.2 mm S1BN1M35 **	 Ø 5 mm S1B1M50 **







Implant systems
Internal hex prosthetic solutions
External hex prosthetic solutions
Conical connection prosthetic solutions
Internal octagon prosthetic solutions
Shapemini prosthetic solutions
MUA and OnePiece components
Accessories, rings and retaining caps
Digital prosthetic solutions
Drills, countersinks taps, stops
Surgical kit
Surgical protocol
Package leaflets



	Internal hex	description   €	
		iMAX NHSI 3.3	iMAX NHSI ShapeOne B Volution SVB
Implant systems			
Internal hex prosthetic solutions		<b>DIRECT IMPRESSION TRANSFER</b> <span>*S1BDTRS screw included   Material Ti-6Al-4V</span>	
External hex prosthetic solutions		Ø 3.5 mm 35.00 €  S1BNDT35 *	 S1BDT35 *
Conical connection prosthetic solutions		Ø 5 mm 35.00 €	 S1BDT50 *
Internal octagon prosthetic solutions		Multifunction abutment 55.00 € Ø 3.5 mm S1BN1A35L *	
Shapemini prosthetic solutions		Multifunction abutment 55.00 € Ø 3.2 mm S1BN1M35L *	
MUA and OnePiece components			
Accessories, rings and retaining caps			
Digital prosthetic solutions		3 components Ø 5 mm 45.00 €  S1BNDT403	 S1BNDT503
Drills, countersinks taps, stops			
Surgical kit			
Surgical protocol			
Package leaflets			

TEMPORARY STRAIGHT ABUTMENTS

S1BRS1 screw included

<div>in PEEK</div> <div>55.00 €</div>	<div>  </div> <div> <div>Ø 3.5 mm</div> <div>S1BN135PP</div> </div>	<div>  </div> <div> <div>Ø 4 mm</div> <div>S1B140PP</div> </div>
<div>Anti rotation in Ti-6Al-4V</div> <div>45.00 €</div>	<div>  </div> <div> <div>Ø 3.5 mm</div> <div>S1BN135P</div> </div>	<div>  </div> <div> <div>Ø 4 mm</div> <div>S1B140P</div> </div>
<div>Rotating in Ti-6Al-4V</div> <div>45.00 €</div>	<div>  </div> <div> <div>Ø 3.5 mm</div> <div>S1BN135PR</div> </div>	<div>  </div> <div> <div>Ø 4 mm</div> <div>S1B140PR</div> </div>












ABUTMENTS FOR WELDED TECHNIQUE

S1BRS1 screw included

<div>55.00 €</div>	<div>  </div> <div>S1BN1TS</div>	<div>  </div> <div>S1B1TS</div>
--------------------	---	--

Implant systems
Internal hex prosthetic solutions
External hex prosthetic solutions
Conical connection prosthetic solutions
Internal octagon prosthetic solutions
Shapemini prosthetic solutions
MUA and OnePiece components
Accessories, rings and retaining caps
Digital prosthetic solutions
Drills, countersinks taps, stops
Surgical kit
Surgical protocol
Package leaflets









Implant systems	Internal hex	iMAX NHSI 3.3		iMAX NHSI ShapeOne B Volution SVB													
		description   €															
Internal hex prosthetic solutions	DEFINITIVE STRAIGHT ABUTMENTS   ANATOMIC																
External hex prosthetic solutions	S1BRS1 screw included   Material Ti-6Al-4V																
Conical connection prosthetic solutions	Friction Fit Ø 3.5 mm  65.00 €	 S1BN135FF	 S1B135F														
Internal octagon prosthetic solutions	Friction Fit Ø 5 mm  65.00 €		 S1B150F														
Shapemini prosthetic solutions	Friction Fit without emergence profile Ø 4.5 mm  65.00 €		 S1B145FS														
MUA and OnePiece components	Friction Fit without emergence profile Ø 5.5 mm  65.00 €		 S1B155FS														
Accessories, rings and retaining caps	Anatomic without Friction Fit  55.00 €	 Ø 4 mm S1BN1 [X] 40	<table><tr><th>Code</th><th>h</th></tr><tr><td>S1BN1140</td><td>1 mm</td></tr><tr><td>S1BN1340</td><td>3 mm</td></tr></table>	Code	h	S1BN1140	1 mm	S1BN1340	3 mm	 Ø 4.5 mm S1B1 [X] 45	<table><tr><th>Code</th><th>h</th></tr><tr><td>S1B1145</td><td>1 mm</td></tr><tr><td>S1B1345</td><td>3 mm</td></tr></table>	Code	h	S1B1145	1 mm	S1B1345	3 mm
Code	h																
S1BN1140	1 mm																
S1BN1340	3 mm																
Code	h																
S1B1145	1 mm																
S1B1345	3 mm																
Digital prosthetic solutions	Anatomic with Friction Fit  65.00 €	 Ø 4 mm S1BN1 [X] 40F	<table><tr><th>Code</th><th>h</th></tr><tr><td>S1BN1140F</td><td>1 mm</td></tr><tr><td>S1BN1340F</td><td>3 mm</td></tr></table>	Code	h	S1BN1140F	1 mm	S1BN1340F	3 mm	 Ø 4.5 mm S1B1 [X] 45F	<table><tr><th>Code</th><th>h</th></tr><tr><td>S1B1145F</td><td>1 mm</td></tr><tr><td>S1B1345F</td><td>3 mm</td></tr></table>	Code	h	S1B1145F	1 mm	S1B1345F	3 mm
Code	h																
S1BN1140F	1 mm																
S1BN1340F	3 mm																
Code	h																
S1B1145F	1 mm																
S1B1345F	3 mm																
Drills, countersinks taps, stops		STICKING BASES															
Surgical kit		S1BRS1 screw included   Material Ti-6Al-4V															
Package leaflets	Friction Fit with emergence profile  65.00 €	 Ø 3.5 mm S1BN135F	 Ø 4.5 mm S1B140F														

24



STICKING BASES *(continued)*

S1BRS1 screw included | Material TI-6Al-4V

Rotating with emergence profile 55.00 €	 Ø 3.5 mm S1BN135R	 Ø 4.5 mm S1B140R
Friction Fit without emergence profile 65.00 €	 Ø 3.5 mm S1BN135FS	 Ø 4.5 mm S1B140FS
Rotating without emergence profile 55.00 €	 Ø 3.5 mm S1BN135RS	 Ø 4.5 mm S1B140RS

ANGLED ABUTMENTS

S1BRS1 screw included | Material TI-6Al-4V


Friction Fit 20° Ø 4.5 mm 75.00 €		S1B250F
Friction Fit without emergence profile 20° Ø 4.5 mm 75.00 €		S1B245FS

ANATOMIC ANGLED ABUTMENTS

S1BRS1 screw included | Material TI-6Al-4V

Without Friction Fit 15°


65.00 €



Ø 4 mm

S1BN2 [X] 1540

Code	h
S1BN211540	1 mm
S1BN231540	3 mm



Ø 4.5 mm







S1B2 [X] 1545

Code	h
S1B211545	1 mm
S1B231545	3 mm









### ANATOMIC ANGLED ABUTMENTS (continued)

S1BRS1 screw included | Material Ti-6Al-4V

<p>Without Friction Fit 25°</p> <p>65.00 €</p>	<div data-bbox="454 390 547 502">  </div> <div data-bbox="479 516 553 541"> <p>Ø 4 mm</p> </div> <div data-bbox="443 548 591 574"> <p><b>S1BN2 [X] 2540</b></p> </div> <table data-bbox="695 479 886 569"> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1BN212540</td><td>1 mm</td></tr> <tr> <td>S1BN232540</td><td>3 mm</td></tr> </table>	Code	h	S1BN212540	1 mm	S1BN232540	3 mm	<div data-bbox="1003 390 1096 502">  </div> <div data-bbox="1016 516 1109 541"> <p>Ø 4.5 mm</p> </div> <div data-bbox="998 548 1127 574"> <p><b>S1B2 [X] 2545</b></p> </div> <table data-bbox="1239 479 1430 569"> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1B212545</td><td>1 mm</td></tr> <tr> <td>S1B232545</td><td>3 mm</td></tr> </table>	Code	h	S1B212545	1 mm	S1B232545	3 mm
Code	h													
S1BN212540	1 mm													
S1BN232540	3 mm													
Code	h													
S1B212545	1 mm													
S1B232545	3 mm													
<p>Friction Fit 15°</p> <p>75.00 €</p>	<div data-bbox="461 652 540 759">  </div> <div data-bbox="479 778 553 803"> <p>Ø 4 mm</p> </div> <div data-bbox="436 810 597 835"> <p><b>S1BN2 [X] 1540F</b></p> </div> <table data-bbox="695 739 886 831"> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1BN211540F</td><td>1 mm</td></tr> <tr> <td>S1BN231540F</td><td>3 mm</td></tr> </table>	Code	h	S1BN211540F	1 mm	S1BN231540F	3 mm	<div data-bbox="1003 652 1084 759">  </div> <div data-bbox="1016 778 1109 803"> <p>Ø 4.5 mm</p> </div> <div data-bbox="992 810 1133 835"> <p><b>S1B2 [X] 1545F</b></p> </div> <table data-bbox="1239 739 1430 831"> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1B211545F</td><td>1 mm</td></tr> <tr> <td>S1B231545F</td><td>3 mm</td></tr> </table>	Code	h	S1B211545F	1 mm	S1B231545F	3 mm
Code	h													
S1BN211540F	1 mm													
S1BN231540F	3 mm													
Code	h													
S1B211545F	1 mm													
S1B231545F	3 mm													
<p>Friction Fit 25°</p> <p>75.00 €</p>	<div data-bbox="454 911 547 1019">  </div> <div data-bbox="479 1035 553 1060"> <p>Ø 4 mm</p> </div> <div data-bbox="436 1067 597 1092"> <p><b>S1BN2 [X] 2540F</b></p> </div> <table data-bbox="695 1001 886 1092"> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1BN212540F</td><td>1 mm</td></tr> <tr> <td>S1BN232540F</td><td>3 mm</td></tr> </table>	Code	h	S1BN212540F	1 mm	S1BN232540F	3 mm	<div data-bbox="1003 911 1096 1019">  </div> <div data-bbox="1016 1035 1109 1060"> <p>Ø 4.5 mm</p> </div> <div data-bbox="992 1067 1133 1092"> <p><b>S1B2 [X] 2545F</b></p> </div> <table data-bbox="1239 1001 1430 1092"> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1B212545F</td><td>1 mm</td></tr> <tr> <td>S1B232545F</td><td>3 mm</td></tr> </table>	Code	h	S1B212545F	1 mm	S1B232545F	3 mm
Code	h													
S1BN212540F	1 mm													
S1BN232540F	3 mm													
Code	h													
S1B212545F	1 mm													
S1B232545F	3 mm													

## CASTABLE ABUTMENTS

S1BRS1 screw included

<p>Gold base POM-C / AU&amp;PGM</p> <p>110.00 €</p>	 <p>Ø 4.5 mm</p> <p>S1BN3GA35 *</p>	 <p>Ø 5 mm</p> <p>S1B3GA50 *</p>
<p>*Gold alloy AU60% PD15% PT24,9% IR0,1% AU e PGM100%. Melting range C° 1350/1460 Tensile modulus GPa 110. Elastic limit MPa 450-720. Stretching %18-12. Breaking load MPa 580-810. Vickers hardness HV5/30 105-205-230.</p>		
<p>POM-C</p> <p>29.00 €</p>	 <p>Ø 3.5 mm</p> <p>S1BN3PC35</p>	 <p>Ø 4.5 mm</p> <p>S1B3PC45</p>
<p>Rotating POM-C</p> <p>29.00 €</p>	 <p>Ø 3.5 mm</p> <p>S1BN3PCR35</p>	 <p>Ø 4.5 mm</p> <p>S1B3PCR45</p>

Internal  
hex

iMAX NHSI 3.3

iMAX NHSI  
ShapeOne B  
Volution SVB

description | €

### CASTABLE ABUTMENTS *(continued)*

S1BRS1 screw included

Titanium base  
Ti-6Al-4V

75.00 €



S1BN3PTC45



S1B3PTC45

### STRAIGHT MUA ABUTMENTS

Mounter included | TIN Treatment on the gold part | Material Ti-6Al-4V

Ø 5 mm

55.00 €



S1BN4 [X]

Code	h
S1BN41	1 mm
S1BN42	2 mm
S1BN43	3 mm
S1BN44	4 mm
S1BN45	5 mm
S1BN46	6 mm



S1B4 [X]

Code	h
S1B41	1 mm
S1B42	2 mm
S1B43	3 mm
S1B44	4 mm
S1B45	5 mm
S1B46	6 mm

### ANGLED MUA ABUTMENTS

Mounter and S1BRS2 screw included | TIN Treatment on the gold part | Material Ti-6Al-4V

18°  
Ø 5 mm  
h 0/2 mm

85.00 €



S1BN518



S1B518

30°  
Ø 5 mm

85.00 €



S1BN53 [X]

Code	h
S1BN532	0/2 mm
S1BN534	2/4 mm



S1B53 [X]

Code	h
S1B532	0/2 mm
S1B534	2/4 mm

### BALL ABUTMENTS

Cah and Calt included | TIN Treatment on the gold part | Material Ti-6Al-4V

Ø 4 mm

45.00 €



S1BN6 [X]

Code	h
S1BN61	1 mm
S1BN62	2 mm
S1BN63	3 mm
S1BN64	4 mm
S1BN65	5 mm



S1B6 [X]

Code	h
S1B61	1 mm
S1B62	2 mm
S1B63	3 mm
S1B64	4 mm
S1B65	5 mm

### IRETOR

TIN Treatment on the gold part

TIN

95.00 €



S1BN8 [XX]

Code	h
S1BN80	0 mm
S1BN81	1 mm
S1BN825	2.5 mm
S1BN835	3.5 mm
S1BN845	4.5 mm
S1BN865	6.5 mm



S1B8 [XX]

Code	h
S1B80	0 mm
S1B81	1 mm
S1B825	2.5 mm
S1B835	3.5 mm
S1B845	4.5 mm
S1B865	6.5 mm

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps








Digital prosthetic  
solutions

Drills, countersinks  
taps, stops

Surgical  
kit

Surgical  
protocol

Package  
leaflets



Implant systems	External hex	iMAX NHSE 3.3		iMAX NHSE									
		description   €											
Internal hex prosthetic solutions	External hex prosthetic solutions	SURGICAL SCREWS											
		Provided with the implant, available as spare parts   Material Ti-6Al-4V											
		Ø 3.5 mm thread 1.6 mm  15.00 €	  S1EHNCS										
		Ø 4.1 mm thread 2 mm  15.00 €		  S1EHCS									
	Internal octagon prosthetic solutions	PROSTHETIC SCREWS											
		Material Ti-6Al-4V											
Shapemini prosthetic solutions	MUA and OnePiece components	For abutments/ lab work thread 1.6 mm  15.00 €	  S1EHNRS1										
		For abutments/ lab work thread 2 mm  15.00 €		  S1EHR1									
	Accessories, rings and retaining caps	Long for transfer thread 1.6 mm  15.00 €	  S1EHNDTRS										
		Long for transfer thread 2 mm  15.00 €		  S1EHDTRS									
	Surgical kit	HEALING SCREWS											
		Material Ti-6Al-4V											
Surgical protocol	Package leaflets	Ø 3.5 mm  25.00 €	  S1EHN35 [x] HC	<table><tr><th>Code</th><th>h</th></tr><tr><td>S1EHN3530HC</td><td>3 mm</td></tr><tr><td>S1EHN3545HC</td><td>4.5 mm</td></tr><tr><td>S1EHN3560HC</td><td>6 mm</td></tr></table>	Code	h	S1EHN3530HC	3 mm	S1EHN3545HC	4.5 mm	S1EHN3560HC	6 mm	
		Code	h										
S1EHN3530HC	3 mm												
S1EHN3545HC	4.5 mm												
S1EHN3560HC	6 mm												

28







HEALING SCREWS *(continued)*










Material Ti-6Al-4V

<div>Ø 4.1 mm</div> <div>25.<sup>00</sup> €</div>	<div>  <div>S1EH41 [xx] HC</div> </div> <div> <table> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1EH4130HC</td><td>3 mm</td></tr> <tr> <td>S1EH4145HC</td><td>5 mm</td></tr> <tr> <td>S1EH4160HC</td><td>7 mm</td></tr> </table> </div>	Code	h	S1EH4130HC	3 mm	S1EH4145HC	5 mm	S1EH4160HC	7 mm
Code	h								
S1EH4130HC	3 mm								
S1EH4145HC	5 mm								
S1EH4160HC	7 mm								
<div>Ø 5 mm</div> <div>25.<sup>00</sup> €</div>	<div>  <div>S1EH50 [xx] HC</div> </div> <div> <table> <tr> <th>Code</th><th>h</th></tr> <tr> <td>S1EH5030HC</td><td>3 mm</td></tr> <tr> <td>S1EH5045HC</td><td>5 mm</td></tr> <tr> <td>S1EH5060HC</td><td>7 mm</td></tr> </table> </div>	Code	h	S1EH5030HC	3 mm	S1EH5045HC	5 mm	S1EH5060HC	7 mm
Code	h								
S1EH5030HC	3 mm								
S1EH5045HC	5 mm								
S1EH5060HC	7 mm								

ANALOG

Material Ti-6Al-4V



<div>21.<sup>00</sup> €</div>	<div>  <div>S1EHNIA</div> </div>	<div>  <div>S1EHIA</div> </div>
<div>Short analog for guided surgery</div> <div>21.<sup>00</sup> €</div>	<div>  <div>S1EHNIA S</div> </div>	<div>  <div>S1EHIA S</div> </div>

	External hex	iMAX NHSE 3.3	iMAX NHSE
		description   €	
Implant systems		DIRECT TRANSFER IMPRESSION *S1EHNDTRS screw included**S1EHDTRS screw included   Material Ti-6Al-4V	
Internal hex prosthetic solutions		Ø 3.5 mm 35.00 €  S1EHNDT35 *	
External hex prosthetic solutions			
Conical connection prosthetic solutions		Multifunction Ø 3.5 mm 65.00 €  S1EHN1A35L *	
Internal octagon prosthetic solutions		Ø 4.1 mm 35.00 €	 S1EHDT41 **
Shapemini prosthetic solutions			
MUA and OnePiece components		Ø 5 mm 35.00 €	 S1EHDT50 **
Accessories, rings and retaining caps		Multifunction 55.00 €	 S1EH1A50L **
Digital prosthetic solutions			
Drills, countersinks taps, stops		TEMPORARY STRAIGHT ABUTMENTS *S1EHNRS1 screw included **S1EHRS1 screw included   Material Ti-6Al-4V	
Surgical kit		Anti rotation 45.00 €  Ø 3.5 mm S1EHN135P *	 Ø 4.1 mm S1EH141P **
Surgical protocol			
Package leaflets		Rotating 45.00 €  Ø 3.5 mm S1EHN135PR *	 Ø 4.1 mm S1EH141PR **

30







ABUTMENTS FOR WELDED TECHNIQUE

\*S1EHNRS1 screw included \*\*S1EHR51 screw included | Material TI-6Al-4V

colored internal prosthetic part	 S1EHN1TS *	 S1EH1TS **
-------------------------------------	---	---





DEFINITIVE STRAIGHT ABUTMENTS

\*S1EHNRS1 screw included \*\*S1EHR51 screw included | Material TI-6Al-4V

Multifunction <i>colored internal prosthetic part</i>	 Ø 3.5 mm S1EHN1A35 *	 Ø 5 mm S1EH1A50 **												
Anti rotation <i>colored internal prosthetic part</i>	 Ø 4 mm S1EHN140 *	 Ø 5 mm S1EH150 **												
Anatomic <i>colored internal prosthetic part</i>	 h Ø 3.5 mm S1EHN1 [X] 35 * <table><tr><th>Code</th><th>h</th></tr><tr><td>S1EHN1135</td><td>1 mm</td></tr><tr><td>S1EHN1335</td><td>3 mm</td></tr></table>	Code	h	S1EHN1135	1 mm	S1EHN1335	3 mm	 h Ø 4.1 mm S1EH1 [X] 41 ** <table><tr><th>Code</th><th>h</th></tr><tr><td>S1EH1141</td><td>1 mm</td></tr><tr><td>S1EH1341</td><td>3 mm</td></tr></table>	Code	h	S1EH1141	1 mm	S1EH1341	3 mm
Code	h													
S1EHN1135	1 mm													
S1EHN1335	3 mm													
Code	h													
S1EH1141	1 mm													
S1EH1341	3 mm													

STICKING BASES

\*S1EHNRS1 screw included \*\*S1EHR51 screw included | Material TI-6Al-4V

Anti rotation colored internal prosthetic part	 Ø 4 mm S1EHN135 *	 Ø 4.1 mm S1EH141 **
Rotating colored internal prosthetic part	 Ø 4 mm S1EHN135R *	 Ø 4.1 mm S1EH141R **




Implant systems	External hex	iMAX NHSE 3.3		iMAX NHSE												
		description   €														
Internal hex prosthetic solutions	External hex prosthetic solutions	ANGLED ABUTMENTS														
		*S1EHNRS1 screw included **S1EHR51 screw included   Material TI-6Al-4V														
Conical connection prosthetic solutions	Internal octagon prosthetic solutions	15° Angled <i>colored internal prosthetic part</i>	 <table><tr><th>Code</th><th>h</th></tr><tr><td>S1EHN211540</td><td>1 mm</td></tr><tr><td>S1EHN231540</td><td>3 mm</td></tr></table> <p>Ø 4 mm S1EHN2 [X] 1540 *</p>	Code	h	S1EHN211540	1 mm	S1EHN231540	3 mm	 <table><tr><th>Code</th><th>h</th></tr><tr><td>S1EH211550</td><td>1 mm</td></tr><tr><td>S1EH231550</td><td>3 mm</td></tr></table> <p>Ø 5 mm S1EH2 [X] 1550 **</p>	Code	h	S1EH211550	1 mm	S1EH231550	3 mm
		Code	h													
S1EHN211540	1 mm															
S1EHN231540	3 mm															
Code	h															
S1EH211550	1 mm															
S1EH231550	3 mm															
Shapenini prosthetic solutions	MUA and OnePiece components	25° Angled <i>colored internal prosthetic part</i>	 <table><tr><th>Code</th><th>h</th></tr><tr><td>S1EH212540</td><td>1 mm</td></tr><tr><td>S1EH232540</td><td>3 mm</td></tr></table> <p>Ø 4 mm S1EHN2 [X] 2540 *</p>	Code	h	S1EH212540	1 mm	S1EH232540	3 mm	 <table><tr><th>Code</th><th>h</th></tr><tr><td>S1EH212550</td><td>1 mm</td></tr><tr><td>S1EH232550</td><td>3 mm</td></tr></table> <p>Ø 5 mm S1EH2 [X] 2550 **</p>	Code	h	S1EH212550	1 mm	S1EH232550	3 mm
		Code	h													
S1EH212540	1 mm															
S1EH232540	3 mm															
Code	h															
S1EH212550	1 mm															
S1EH232550	3 mm															
Accessories, rings and retaining caps	Digital prosthetic solutions	CASTABLE ABUTMENTS														
		*S1EHNRS1 screw included **S1EHR51 screw included   Material POM-C														
Drills, countersinks taps, stops	Surgical kit	Gold base POM-C / AU&PGM	 <p>Ø 3.5 mm S1EHN3GA *</p>	 <p>Ø 4.1 mm S1EH3GA **</p>												
		Gold alloy AU60% PD15% PT24,9% IR0,1% AU e PGM100%. Melting range C° 1350/1460 Tensile modulus GPa 110. Elastic limit MPa 450-720. Stretching %18-12. Breaking load MPa 580-810. Vickers hardness HV5/30 105-205-230.														
Surgical protocol	Package leaflets	Gold base Rotating POM-C / AU&PGM	 <p>Ø 3.5 mm S1EHN3GAR *</p>	 <p>Ø 4.1 mm S1EH3GAR **</p>												
		Gold alloy AU60% PD15% PT24,9% IR0,1% AU e PGM100%. Melting range C° 1350/1460 Tensile modulus GPa 110. Elastic limit MPa 450-720. Stretching %18-12. Breaking load MPa 580-810. Vickers hardness HV5/30 105-205-230.														
			 <p>Ø 3.5 mm S1EHN3PC35 *</p>	 <p>Ø 4.1 mm S1EH3PC41 *</p>												
		Rotating	 <p>Ø 3.5 mm S1EHN3PCR35 *</p>	 <p>Ø 4.1 mm S1EH3PCR41 *</p>												

32

STRAIGHT MUA ABUTMENTS

Mounter included | TIN Treatment on the gold part | Material TI-6Al-4V


55.00 €



Ø 3.5 mm

S1EHN4 [X]

Code	h
S1EHN41	1 mm
S1EHN42	2 mm
S1EHN43	3 mm
S1EHN44	4 mm
S1EHN45	5 mm







Ø 4.1 mm

S1EH4 [X]

Code	h
S1EH41	1 mm
S1EH42	2 mm
S1EH43	3 mm
S1EH44	4 mm
S1EH45	5 mm

ANGLED MUA ABUTMENTS


\*S1EHNRS2 screw included \*\*S1EHS2 screw included | Mounter included  
TIN Treatment on the gold part | Material TI-6Al-4V

18° h 0/2 mm  85.00 €	 Ø 3.5 mm S1EHN518 *	 Ø 4.1 mm S1EH518 **						
30°  85.00 €	 Ø 3.5 mm - h 0/2 mm S1EHN532 *	 Ø 4.1 mm S1EH53[X]** <table border="1" data-bbox="1206 1157 1391 1242"><thead><tr><th>Code</th><th>h</th></tr></thead><tbody><tr><td>S1EH532</td><td>0/2 mm</td></tr><tr><td>S1EH534</td><td>2/4 mm</td></tr></tbody></table>	Code	h	S1EH532	0/2 mm	S1EH534	2/4 mm
Code	h							
S1EH532	0/2 mm							
S1EH534	2/4 mm							

BALL ABUTMENTS


Cah and Calt included | TIN Treatment on the gold part | Material TI-6Al-4V

45.00 €



Ø 3.5 mm  
S1EHN6 [X]

Code	h
S1EHN61	1 mm
S1EHN62	2 mm
S1EHN63	3 mm
S1EHN64	4 mm



Ø 4.1 mm  
S1EH6 [X]


Code	h
S1EH61	1 mm
S1EH62	2 mm
S1EH63	3 mm
S1EH64	4 mm

iRETOR

TIN Treatment on the gold part


TIN

95.00 €





S1EHN8 [XX]

Code	h
S1EHN80	0 mm
S1EHN81	1 mm
S1EEHN825	2.5 mm
S1EHN835	3.5 mm
S1EHN845	4.5 mm
S1EHN865	6.5 mm



S1EH8 [XX]

Code	h
S1EH80	0 mm
S1EH81	1 mm
S1EH825	2.5 mm
S1EH835	3.5 mm
S1EH845	4.5 mm
S1EH865	6.5 mm

Implant systems	Conical connection	iMAX NHSIC Narrow		iMAX NHSIC Regular	
		description   €			
Internal hex prosthetic solutions	External hex prosthetic solutions	SURGICAL SCREWS			
		Provided with the implant, available as spare parts   Material Ti-6Al-4V			
		15.00 €			
			NHSICNCS		NHSICCS
Conical connection prosthetic solutions		PROSTHETIC SCREWS			
		Material Ti-6Al-4V			
Internal octagon prosthetic solutions	Shapemini prosthetic solutions	For abutments/ lab work thread 1.72 mm			
		15.00 €		S1BRS1	
MUA and OnePiece components	Accessories, rings and retaining caps	For comp MUA /lab work thread 1.72 mm			
		15.00 €		S1BRS2	
Digital prosthetic solutions		Long for transfer thread 1.72 mm			
		15.00 €		S1BDTRS	
Drills, countersinks taps, stops	Surgical kit	Long for transfer MUA thread 1.4 mm			
		16.00 €		S1BDTRSA	
Surgical protocol	Package leaflets	ANALOG			
		Material Ti-6Al-4V			
		26.00 €			
			Ø 4 mm NHSICNIA		Ø 4.5 mm NHSICIA
		Guided surgery			
		30.00 €	Ø 4 mm NHSICNIAS		Ø 4.5 mm NHSICIAS


34



HEALING SCREWS

S1BRS1 screw included | Material TI-6Al-4V

30.00 €




h

Ø 4 mm

NHSICN40 [x] HC

Code	h
NHSICN4030HC	3 mm
NHSICN4045HC	4.5 mm
NHSICN4060HC	6 mm



h

Ø 4.5 mm


NHSIC45 [x] HC

Code	h
NHSIC4530HC	3 mm
NHSIC4545HC	4.5 mm
NHSIC4560HC	6 mm

DIRECT TRANSFER IMPRESSION

S1BDTRS screw included | Material TI-6Al-4V

45.00 €




h

Ø 4 mm

NHSICNDT [x] 40

Code	h
NHSICNDT140	1.5 mm
NHSICNDT240	2.5 mm
NHSICNDT340	4 mm



h

Ø 4.5 mm


NHSICDT [x] 45

Code	h
NHSICDT145	1.5 mm
NHSICDT245	2.5 mm
NHSICDT345	4 mm

PEEK TEMPORARY ABUTMENTS

S1BRS1 screw included | Material TI-6Al-4V

50.00 €




h

Ø 4 mm

NHSICN1 [x] 40PP

Code	h
NHSICN1140PP	1.5 mm
NHSICN1240PP	2.5 mm
NHSICN1440PP	4 mm



h

Ø 4.5 mm


NHSIC1 [x] 45PP

Code	h
NHSIC1145PP	1.5 mm
NHSIC1245PP	2.5 mm
NHSIC1445PP	4 mm

STRAIGHT ABUTMENTS

S1BRS1 screw included | Material TI-6Al-4V

85.00 €




h

Ø 4 mm

NHSICN1 [x] 40

Code	h
NHSICN1140	1.5 mm
NHSICN1240	2.5 mm
NHSICN1440	4 mm













h

Ø 4.5 mm








NHSIC1 [x] 45

Code	h
NHSIC1145	1.5 mm
NHSIC1245	2.5 mm
NHSIC1445	4 mm

Implant systems	Conical connection	iMAX NHSIC Narrow		iMAX NHSIC Regular																
		description   €																		
Internal hex prosthetic solutions	External hex prosthetic solutions	STICKING BASES																		
		S1BRS1 screw included   Material Ti-6Al-4V																		
Conical connection prosthetic solutions	Internal octagon prosthetic solutions	75.00 €	 <p>Ø 4 mm NHSICN1 [X] 40SB</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN1140SB</td><td>1.5 mm</td></tr><tr><td>NHSICN1240SB</td><td>2.5 mm</td></tr><tr><td>NHSICN1440SB</td><td>4 mm</td></tr></table>	Code	h	NHSICN1140SB	1.5 mm	NHSICN1240SB	2.5 mm	NHSICN1440SB	4 mm	 <p>Ø 4.5 mm NHSIC1 [X] 45SB</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSIC1145SB</td><td>1.5 mm</td></tr><tr><td>NHSIC1245SB</td><td>2.5 mm</td></tr><tr><td>NHSIC1445SB</td><td>4 mm</td></tr></table>	Code	h	NHSIC1145SB	1.5 mm	NHSIC1245SB	2.5 mm	NHSIC1445SB	4 mm
		Code	h																	
NHSICN1140SB	1.5 mm																			
NHSICN1240SB	2.5 mm																			
NHSICN1440SB	4 mm																			
Code	h																			
NHSIC1145SB	1.5 mm																			
NHSIC1245SB	2.5 mm																			
NHSIC1445SB	4 mm																			
Shapemini prosthetic solutions	MUA and OnePiece components	Rotating 75.00 €	 <p>Ø 4 mm NHSICN1 [X] 40RSB</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN1140RSB</td><td>1.5 mm</td></tr><tr><td>NHSICN1240RSB</td><td>2.5 mm</td></tr><tr><td>NHSICN1440RSB</td><td>4 mm</td></tr></table>	Code	h	NHSICN1140RSB	1.5 mm	NHSICN1240RSB	2.5 mm	NHSICN1440RSB	4 mm	 <p>Ø 4.5 mm NHSIC1 [X] 45RSB</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSIC1145RSB</td><td>1.5 mm</td></tr><tr><td>NHSIC1245RSB</td><td>2.5 mm</td></tr><tr><td>NHSIC1445RSB</td><td>4 mm</td></tr></table>	Code	h	NHSIC1145RSB	1.5 mm	NHSIC1245RSB	2.5 mm	NHSIC1445RSB	4 mm
		Code	h																	
NHSICN1140RSB	1.5 mm																			
NHSICN1240RSB	2.5 mm																			
NHSICN1440RSB	4 mm																			
Code	h																			
NHSIC1145RSB	1.5 mm																			
NHSIC1245RSB	2.5 mm																			
NHSIC1445RSB	4 mm																			
Digital prosthetic solutions	Drills, countersinks taps, stops	ANGLED ABUTMENTS																		
		S1BRS1 screw included   Material Ti-6Al-4V																		
Surgical kit	Surgical protocol	15° Angled 85.00 €	 <p>Ø 4 mm NHSICN2 [X] 1540</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN211540</td><td>1.5 mm</td></tr><tr><td>NHSICN221540</td><td>2.5 mm</td></tr><tr><td>NHSICN241540</td><td>4 mm</td></tr></table>	Code	h	NHSICN211540	1.5 mm	NHSICN221540	2.5 mm	NHSICN241540	4 mm	 <p>Ø 4.5 mm NHSIC2 [X] 1545</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSIC211545</td><td>1.5 mm</td></tr><tr><td>NHSIC221545</td><td>2.5 mm</td></tr><tr><td>NHSIC241545</td><td>4 mm</td></tr></table>	Code	h	NHSIC211545	1.5 mm	NHSIC221545	2.5 mm	NHSIC241545	4 mm
		Code	h																	
NHSICN211540	1.5 mm																			
NHSICN221540	2.5 mm																			
NHSICN241540	4 mm																			
Code	h																			
NHSIC211545	1.5 mm																			
NHSIC221545	2.5 mm																			
NHSIC241545	4 mm																			
Package leaflets	POM-C/ AU&PGM	25° Angled 85.00 €	 <p>Ø 4 mm NHSICN2 [X] 2540</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN212540</td><td>1.5 mm</td></tr><tr><td>NHSICN222540</td><td>2.5 mm</td></tr><tr><td>NHSICN242540</td><td>4 mm</td></tr></table>	Code	h	NHSICN212540	1.5 mm	NHSICN222540	2.5 mm	NHSICN242540	4 mm	 <p>Ø 4.5 mm NHSIC2 [X] 2545</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSIC212545</td><td>1.5 mm</td></tr><tr><td>NHSIC222545</td><td>2.5 mm</td></tr><tr><td>NHSIC242545</td><td>4 mm</td></tr></table>	Code	h	NHSIC212545	1.5 mm	NHSIC222545	2.5 mm	NHSIC242545	4 mm
		Code	h																	
NHSICN212540	1.5 mm																			
NHSICN222540	2.5 mm																			
NHSICN242540	4 mm																			
Code	h																			
NHSIC212545	1.5 mm																			
NHSIC222545	2.5 mm																			
NHSIC242545	4 mm																			
		CASTABLE GOLD BASE ABUTMENTS																		
		S1BRS1 screw included																		
		120.00 €	 <p>Ø 4 mm NHSICN3GA [X] 40</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN3GA140</td><td>1.5 mm</td></tr><tr><td>NHSICN3GA240</td><td>2.5 mm</td></tr><tr><td>NHSICN3GA440</td><td>4 mm</td></tr></table>	Code	h	NHSICN3GA140	1.5 mm	NHSICN3GA240	2.5 mm	NHSICN3GA440	4 mm	 <p>Ø 4.5 mm NHSIC3GA [X] 45</p> <table><tr><th>Code</th><th>h</th></tr><tr><td>NHSIC3GA145</td><td>1.5 mm</td></tr><tr><td>NHSIC3GA245</td><td>2.5 mm</td></tr><tr><td>NHSIC3GA445</td><td>4 mm</td></tr></table>	Code	h	NHSIC3GA145	1.5 mm	NHSIC3GA245	2.5 mm	NHSIC3GA445	4 mm
Code	h																			
NHSICN3GA140	1.5 mm																			
NHSICN3GA240	2.5 mm																			
NHSICN3GA440	4 mm																			
Code	h																			
NHSIC3GA145	1.5 mm																			
NHSIC3GA245	2.5 mm																			
NHSIC3GA445	4 mm																			
Gold alloy AU60% PD15% PT24,9% IR0,1% AU e PGM100%. Melting range C° 1350/1460 Tensile modulus GPa 110. Elastic limit MPa 450-720. Stretching %18-12. Breaking load MPa 580-810. Vickers hardness HV5/30 105-205-230.																				

36







Implant systems	Internal octagon	ShapeOne T		ShapeOne Tn											
		description  €													
Internal hex prosthetic solutions	External hex prosthetic solutions	SURGICAL SCREWS													
		*Provided with the implant, available as spare parts   S1BRS1 screw included													
Conical connection prosthetic solutions	Internal octagon prosthetic solutions	Surgical screw Ø 4.8 mm thread 2 mm  15.00 €	 S1TCS												
Shapemini prosthetic solutions	Internal octagon prosthetic solutions	PROSTHETIC SCREWS													
		Material TI-6Al-4V													
MU and OnePiece components	Accessories, rings and retaining caps	For abutments/ lab work thread 2 mm  15.00 €	 S1TRS1	 S1TRS2											
Digital prosthetic solutions	Drills, countersinks taps, stops	ANALOG													
		Material TI-6Al-4V													
Surgical kit	Surgical protocol	21.00 €	 S1TIA												
		22.00 €		 S1TNIA [xx]	<table><tr><th>Code</th><th>Ø</th></tr><tr><td>S1TNIA37</td><td>3.7 mm</td></tr><tr><td>S1TNIA41</td><td>4.1 mm</td></tr><tr><td>S1TNIA47</td><td>4.7 mm</td></tr></table>	Code	Ø	S1TNIA37	3.7 mm	S1TNIA41	4.1 mm	S1TNIA47	4.7 mm		
Code	Ø														
S1TNIA37	3.7 mm														
S1TNIA41	4.1 mm														
S1TNIA47	4.7 mm														
Package leaflets		DIRECT IMPRESSION TRANSFER													
		S1BDTRS screw included   Material TI-6Al-4V													
		Mounter transfer straight definitive abutment Ø 3.5 mm  55.00 €	 S1TN1A35L	<table><tr><th>Code</th><th>Ø</th></tr><tr><td>S1TN1A35L</td><td>3.5 mm</td></tr><tr><td>S1TN1A37L</td><td>3.7 mm</td></tr><tr><td>S1TN1A41L</td><td>4.1 mm</td></tr><tr><td>S1TN1A47L</td><td>4.7 mm</td></tr></table>	Code	Ø	S1TN1A35L	3.5 mm	S1TN1A37L	3.7 mm	S1TN1A41L	4.1 mm	S1TN1A47L	4.7 mm	
		Code	Ø												
S1TN1A35L	3.5 mm														
S1TN1A37L	3.7 mm														
S1TN1A41L	4.1 mm														
S1TN1A47L	4.7 mm														
		INDIRECT IMPRESSION TRANSFER													
		S1BRS1 screw included   Material TI-6Al-4V													
		Mounter transfer straight definitive abutment  55.00 €	 S1TN1A [xx]	<table><tr><th>Code</th><th>Ø</th></tr><tr><td>S1TN1A35</td><td>3.5 mm</td></tr><tr><td>S1TN1A37</td><td>3.7 mm</td></tr><tr><td>S1TN1A41</td><td>4.1 mm</td></tr><tr><td>S1TN1A47</td><td>4.7 mm</td></tr></table>	Code	Ø	S1TN1A35	3.5 mm	S1TN1A37	3.7 mm	S1TN1A41	4.1 mm	S1TN1A47	4.7 mm	
Code	Ø														
S1TN1A35	3.5 mm														
S1TN1A37	3.7 mm														
S1TN1A41	4.1 mm														
S1TN1A47	4.7 mm														

38









	Shape Mini	ShapeMini with ball connection in TIN	ShapeMini with fixed abutment
Implant systems	description   €		
Internal hex prosthetic solutions	Material Ti-6Al-4V		
External hex prosthetic solutions	Analog 21.00 €	 SMIA	 with fixed abutment SMIAF
Conical connection prosthetic solutions	Castable POM-C abutment 29.00 €	 SMC	
Internal octagon prosthetic solutions			
Shapemini prosthetic solutions	Straight abutment 29.00 €		 SMF
MUA and OnePiece components			
Accessories, rings and retaining caps			
Digital prosthetic solutions			
Drills, countersinks taps, stops			
Surgical kit			
Surgical protocol			
Package leaflets			

MUA	Onepiece iMAXMUA		MUA abutments	
	description   €			

RETAINING SCREWS

Material TI-6Al-4V

Retaining screw for abutment	 NHSMRS1	 S1BRS3
15.00 €		
Long retaining screw for MUA transfer thread 1.72 mm	 NHSMTRS	 S1BDTRSA
15.00 €		


HEALING CAPS

\*NHSMRS1 screw included | \*\*S1BRS3 screw included | Material TI-6Al-4V

Healing TI-6Al-4V cup	 NHSMHCSRA *	 S1BHCSRA **
25.00 €		
POM-C Healing cup	 NHSMHCSRAA *	 S1BHCSRAA **
25.00 €		

DIRECT IMPRESSION TRANSFER

\*NHSMTRS screw included | \*\*S1BDTRSA screw included | Material TI-6Al-4V

With long screw	 NHSMDCSRA *	 S1BDTCSRA **
35.00 €		

INDIRECT IMPRESSION TRANSFER

\*NHSMRS1 screw included | \*\*S1BRS3 screw included | Material TI-6Al-4V

	 NHSMITCSRA *	 S1BITCSRA **
55.00 €		

Implant systems
Internal hex prosthetic solutions
External hex prosthetic solutions
Conical connection prosthetic solutions
Internal octagon prosthetic solutions
Shapemini prosthetic solutions
MUA and OnePiece components
Accessories, rings and retaining caps
Digital prosthetic solutions
Drills, countersinks taps, stops
Surgical kit
Surgical protocol
Package leaflets

## MUA abutments

description | €

## ANALOG

Material Ti-6Al-4V

21.<sup>00</sup> €



NHSMIASRA



S1BIASRA

## TEMPORARY STRAIGHT ABUTMENTS

\*NHSMRS1 screw included | \*\*S1BRS3 screw included | Material Ti-6Al-4V

45.00 €



NHSMPTTA\*



S1BPTTA \*\*

## DEFINITIVE STRAIGHT ABUTMENTS

\*NHSMRS1 screw included | \*\*S1BRS3 screw included | Material Ti-6Al-4V

55.<sup>00</sup> €



NHSMTTA\*



S1BTTA \*\*

## THREADED ABUTMENTS FOR BAR

\*NHSMRS1 screw included | \*\*S1BRS3 screw included | Material Ti-6Al-4V

55.<sup>00</sup> €



NHSMTTAB\*



S1BTTAB \*\*

## CASTABLE ABUTMENTS

\*NHSMRS1 screw included | \*\*S1BRS3 screw included | Material POM-C

29.00 €



NHSMPPC \*



S1BPCC \*\*

## ABUTMENTS FOR WELDED TECHNIQUE

\*NHSMRS1 screw included | \*\*S1BR53 screw included | Material TI-6Al-4V

55.<sup>00</sup> €



NHSMTS\*



S1BTS \*\*




Rings and Caps

Accessories, Rings and retaining Caps for pro overdenture connection

description | €


CONTAINMENT RING

Material TI-6Al-4V




Containment ring		
15. <sup>00</sup> €	CAH	

CONTAINMENT CAP

Material Nylon

Nylon containment cap		
7. <sup>50</sup> €	CALT	

iRETOR ACCESSORIES


Analog for iRETOR® female coupling Q.tà: 4 Al		
50. <sup>00</sup> €	8530	
Indirect impression tear pin for iRETOR® Q.tà: 4 Al		
50. <sup>00</sup> €	8505	
Set of cups and rings for parallel implants Q.tà: 1 Ti   Nylon		
27. <sup>50</sup> €	S1B85	

Implant systems
Internal hex prosthetic solutions
External hex prosthetic solutions
Conical connection prosthetic solutions
Internal octagon prosthetic solutions
Shapemini prosthetic solutions
MUA and OnePiece components
Accessories, rings and retaining caps
Digital prosthetic solutions
Drills, countersinks taps, stops
Surgical kit
Surgical protocol
Package leaflets

--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

SCREWS and SCREWDRIVERS











Material TI-6Al-4V

15.00 €	<div>  <p>Vite di ritenzione filetto 1.72 mm S1BRS1</p> </div>			
Torx screw 15.00 €	<div>  <p>internal hex S1BRS1T6</p> </div>	<div>  <p>external hex S1EHR1T6</p> </div>	<div>  <p>external hex S1ENHR1T6</p> </div>	<div>  <p>internal octagon S1TRS1T6</p> </div>
65.00 €	<div>  <p>Torx screwdriver MDLAA</p> </div>			

DIGITAL SCAN ABUTMENTS

Material Peek

	iMAX NHSI 3.3		iMAX NHSI - ShapeOne B - Volution SVB	
With Titanium base and S1BRS1 screw included  55.00 €	<div></div> S1BNSA		<div></div> S1BSA	
	iMAX NHSE 3.3		iMAX NHSE	
*S1EHNRS1 screw included **S1EHR1S1 screw included  55.00 €	<div></div> S1EHNSA *		<div></div> S1EHSA **	
	iMAX NHSIC Narrow		iMAX NHSIC Regular	
Titanium base S1BRS1 screw included  65.00 €	<div></div> NHSICNSA		<div></div> NHSICSA	
	ShapeOne T			
Titanium base S1TRS1 screw included  55.00 €	<div></div> S1TSA			

	Digital	Prosthetic solutions																	
Implant systems	description   €																		
Internal hex prosthetic solutions	DIGITAL SCAN ABUTMENTS (continued) <span>Materiale Peek</span>																		
External hex prosthetic solutions	<div><div>Onepiece iMAXMUA</div><div>Scan abutment (MUA)</div></div>																		
Conical connection prosthetic solutions	<div><div><div>*NHSMRS1 screw included **S1BRS3 screw included 55.00 €</div><div> NHSMSAA *</div></div><div><div> S1BSAA **</div></div></div>																		
Internal octagon prosthetic solutions	DIGITAL T-BASE <span>Material Ti-6Al-4V</span>																		
Shapenini prosthetic solutions	iMAX NHSI 3.3																		
MUA and Onepiece components	<div><div><div>S1BRS1 screw included 65.00 €</div><div> S1BN1DTB</div></div><div><div> Rotating S1BN1DTBR</div></div></div>																		
Accessories, rings and retaining caps	iMAX NHSI - ShapeOne B - Volution SVB																		
Digital prosthetic solutions	<div><div><div>S1BRS1 screw included 65.00 €</div><div> S1B1DTB</div></div><div><div> Rotating S1B1DTBR</div></div></div>																		
Drills, countersinks taps, stops	iMAX NHSE 3.3																		
Surgical kit	<div><div><div>S1EHNRS1 screw included 65.00 €</div><div> S1EHN1DTB</div></div><div><div> Rotating S1EHN1DTBR</div></div></div>																		
Surgical protocol	iMAX NHSE																		
Package leaflets	<div><div><div>S1BRS1 screw included 75.00 €</div><div> NHSICN [x] DTB</div><div><table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN1DTB</td><td>1 mm</td></tr><tr><td>NHSICN2DTB</td><td>2 mm</td></tr><tr><td>NHSICN4DTB</td><td>4 mm</td></tr></table></div></div><div><div> Rotating NHSICN [x] DTBR</div><div><table><tr><th>Code</th><th>h</th></tr><tr><td>NHSICN1DTBR</td><td>1 mm</td></tr><tr><td>NHSICN2DTBR</td><td>2 mm</td></tr><tr><td>NHSICN4DTBR</td><td>4 mm</td></tr></table></div></div></div>			Code	h	NHSICN1DTB	1 mm	NHSICN2DTB	2 mm	NHSICN4DTB	4 mm	Code	h	NHSICN1DTBR	1 mm	NHSICN2DTBR	2 mm	NHSICN4DTBR	4 mm
Code	h																		
NHSICN1DTB	1 mm																		
NHSICN2DTB	2 mm																		
NHSICN4DTB	4 mm																		
Code	h																		
NHSICN1DTBR	1 mm																		
NHSICN2DTBR	2 mm																		
NHSICN4DTBR	4 mm																		

46



Prosthetic solutions

DIGITAL T-BASE (continued)

Material TI-6Al-4V

iMAX NHSIC Regular

S1BRS1  
screw included  
65.00 €



NHSIC [x] DTB

Code	h
NHSIC1DTB	1 mm
NHSIC2DTB	2 mm
NHSIC4DTB	4 mm



Rotante  
NHSIC [x] DTBR

Code	h
NHSIC1DTBR	1 mm
NHSIC2DTBR	2 mm
NHSIC4DTBR	4 mm

ShapeOne T

S1TRS1  
screw included  
65.00 €



S1T1DTB



Rotating  
S1T1DTBR

descrizione | €

Onepiece iMAXMUA

Monconi a spalla (MUA)

\*NHSMRS1  
screw included  
\*\*S1BRS1  
screw included  
65.00 €



Rotating  
NHSMTTADTB \*



Rotating  
S1BTTADTB \*\*

DIGITAL T-BASE FOR CEREC

Material TI-6Al-4V

iMAX NHSI 3.3

For internal  
hex  
65.00 €



S1BN1DTBC



Rotating  
S1BN1DTBCR

iMAX NHSI - ShapeOne B - Volution SVB

For internal  
hex  
65.00 €



S1B1DTBC



Rotating  
S1B1DTBCR

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps










Digital prosthetic  
solutions

Drills, countersinks  
taps, stops

Surgical  
kit

Surgical  
protocol

Package  
leaflets

	Digital	Prosthetic solutions	
Implant systems	description   €		
Internal hex prosthetic solutions	<div>DIGITAL ANALOG</div> <div>Material Peek</div>		
External hex prosthetic solutions		iMAX NHSI 3.3	iMAX NHSI - ShapeOne B - Volution SVB
Conical connection prosthetic solutions	For internal hex 21.00 €	 S1BNDIA	 S1BDIA
Internal octagon prosthetic solutions	For external hex 21.00 €	 S1EHNDIA	 S1EHDIA
Shapenini prosthetic solutions	For conical connection 26.00 €	iMAX NHSIC Narrow  NHSICNDIA	iMAX NHSIC Regular  NHSICDIA
MUA and OnePiece components		ShapeOne T	
Accessories, rings and retaining caps	For internal octagon 21.00 €	 S1TDIA	
Digital prosthetic solutions		Onepiece iMAXMUA	Scan abutment (MUA)
Drills, countersinks taps, stops	Per MUA 21.00 €	 NHSMDIASRA	 S1BDIASRA
Surgical kit			
Surgical protocol			
Package leaflets			

## One surgical kit for all the implant systems

The image shows the RES Surgical Kit, an open plastic case containing various dental implant drivers and tools. The kit is organized into rows and columns, with labels for different components such as 'PRO. DRIVERS', 'IMPLANT DRIVERS', 'SPARES IMPLANT DRIVERS', and 'PARALLEL PIN'. The case is white with a blue interior lining.

- 1- Dismantle all the compound parts.
- 2- Rinse abundantly with cold or lukewarm water for 2-5 minutes.
- 3- Leave the instruments for 10 minutes in an ultrasonic cleaner with a neutral pH enzymatic detergent diluted in water according to the product instructions.
- 4- Wash the instruments with water for 3 minutes.

The guidelines for sterilization are listed below.  
Exceeding these sterilization limits may cause deterioration of the plastic components.


49

## THE DRILLS







For a proper osteotomy and to maintain the integrity of bone quality, the maximum recommended speed is 800rpm with direct input on drill of saline solution to facilitate cooling. All drills are made of medical stainless steel and subjected to hardening heat treatment. Maximum recommended number of use of the devices 30 times.

The helical milling cutters have reference laser markings that identify the depth to reach, until 8mm with a thin white line, from 10 to 13mm with a white

band in which at half is identified the 11.5mm height, and finally a thin white line for the 16mm. This system gives a clear and intuitive glance of the depth level reached by the drill. 4.5 and 6.5mm are not present to avoid confusion in reading the demarcation lines, and being these measures close to the nerve, it is always recommended to use stop by 4.5 and 6.5mm.



33/38 - 38/44 - 44/48 are flathead drills. They can not be used to cut, but only to help you to insert the implant





## THE COUNTERSINKS

The countersinks are used when there is the need to enlarge the initial part of the hole created to adapt this shape to the neck of the implant to be inserted. The maximum recommended speed is 300rpm with direct input on drill of saline solution to facilitate cooling. The countersink should be used in perfect axis with the osteotomy to avoid its ovalization in the coronal part. The countersinks present two laser markings that

identify the depth to be reached on the basis of the bone consistency, at 1.4mm for a “D3” bone, at 2.8mm for both “D2” and “D1” bones. Above the marking at 2.8mm, the countersink continues with a cylindrical geometry that does not compromise the osteotomy although more deeply inserted.



## THE TAPS

In very dense bone (Type I) it is recommended to use a tap with the same system profile to insert. The tap is sharper than the implant and it allows to prepare the implantation site with reduced trauma.


The maximum recommended speed is 30 rpm with direct input on tap of saline solution to facilitate cooling.



## THE STOPS


The length of the stops ranges from 4.5mm to 13mm and they are available for all the implant lengths. 33/38

- 38/44 - 44/48 are flathead drills. They can not be used to cut, but only to help you to insert the implant.



## RATCHET 10/100 NCM and CONNECTORS

In the kit all the connectors have a handpiece attack that may be used both in manual mode and with the ratchet, thanks to the special washer (WH2). The insertion torque for the immediate loading will be in the range from 35 to 50Ncm. For the conventional load the insertion torque should never exceed 70Ncm.



### RATCHET ADAPTOR

code: **WH2**

Ø 8 mm reinforced that adding solidity



### RATCHET WRENCH

code: **TW2**

Mounted on the reverse ratchet, it allows to measure up to 100Ncm2 without breaking the rod through the stop final race



### REVERSE FIXED RATCHET + RATCHET ADAPTOR + RATCHET WRENCH

code: **TWA2**



Reverse to change direction of unscrewing and screwing without having to remove and replace the ratchet..

1. The new ratchet uses connectors with **HANDLE ATTACK**

2. the **DRIVE HEXAGON** confers greater solidity during the tightening phase

3. **POSITION HEXAGON** during the insertion allows to see the internal positioning of the hexagon

4. **HEXAGON CONNECTION** that reaches the stop of the system, has a steel retention ring reinforced and raised both to avoid interferences during the implant insertion and to reduce the wear of the retention



### REVERSE FIXED RATCHET

code: **R2**

It allows to screw and unscrew without having to pull out and turn the adapter



### REVERSE FIXED RATCHET + RATCHET ADAPTOR

code: **RA2**



Housing for 8mm washer to confer greater resistance to higher torque.



Under the 100Ncm is present a safety catch to prevent the leakage of the dragging arm, avoiding its breakage.

SURGICAL KIT COMPOSITION



Code	Description
DE	drill extender
LD	lance drill
D20M	pilot drill Ø 2.0
D2024M	drill Ø 2.0 2.4 mm
D2428M	drill Ø 2.4 2.8 mm
D2833M	drill Ø 2.8 3.3 mm
D3338M	drill Ø 3.3 3.8 mm
D3844M	drill Ø 3.8 4.4 mm
D4448M	drill Ø 4.4 4.8 mm
CSD33	countersink Ø 3.3
CSD37	countersink Ø 3.7
CSD41	countersink Ø 4.1
CSD47	countersink Ø 4.7
CSD52	countersink Ø 5.2
TAPS1B37	tap Ø 3.7
TAPS1B41	tap Ø 4.1
TAPS1B47	tap Ø 4.7
TAPiMAX33	tap Ø 3.3
TAPiMAX37	tap Ø 3.7
TAPiMAX41	tap Ø 4.1

Code	Description
TAPiMAX47	tap Ø 4.7
TAPiMAX52	tap Ø 5.2
DS43341M	stop 4.5
DS63341M	stop 6.5
DS83341M	stop 8
DS103341M	stop 10
DS113341M	stop 11.5
DS133341M	stop 13
TWA2	complete ratchet
THDDS	short driver
THDDL	long driver
PP	parallel pin 0°
PP18	parallel pin 18°
PP30	parallel pin 30°

INTERNAL HEX DRIVERS

Code	Description	Hex
HDH21S	short driver hex 2.1	internal
HDH21L	long driver hex 2.1	internal
HDH25S	short driver hex 2.5	internal
HDH25L	long driver hex 2.5	internal

EXTERNAL OCTAGON DRIVERS

Code	Description	Octagon
HDH31S	short driver oct. 3.1	external
HDH31L	long driver oct. 3.1	external

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps
















Digital prosthetic  
solutions

Drills, countersinks  
taps, stops

Surgical  
kit

Surgical  
protocol

Package  
leaflets

Implant systems	Surgical Kit	Drills, Countersinks, Taps, Stops						
		description   €						
Internal hex prosthetic solutions	External hex prosthetic solutions	DRILL EXTENDER					Material Surgical steel	
		40.00 €						DE
Conical connection prosthetic solutions	Shapemini prosthetic solutions	LANCE DRILL					Material Surgical steel	
Internal octagon prosthetic solutions		Ø 2 mm 45.00 €						LD
MUA and OnePiece components	Accessories, rings and retaining caps	TISSUE PUNCH					Material Surgical steel	
		48.00 €					Ø 3.3 mm HDHCSN	Ø 4 mm HDHCSR
Digital prosthetic solutions	Drills, countersinks taps, stops	DRILLS					Material Surgical steel	
		Medium h 36 mm 70.00 €						
	Surgical Kit	Ø 2 mm D20M	Ø 2.0 2.4 mm D2024M	Ø 2.4 2.8 mm D2428M	Ø 2.8 3.3 mm D2833M	Ø 3.3 3.8 mm D3338M	Ø 3.8 4.4 mm D3844M	Ø 4.4 4.8 mm D4448M
Surgical protocol	Package leaflets	COUNTERSINKS					Material Surgical steel	
		90.00 €						
		Ø 3.3 mm CSD33	Ø 3.7 mm CSD37	Ø 4.1 mm CSD41	Ø 4.7 mm CSD47	Ø 5.2 mm CSD52		

54

Drills, Countersinks, Taps, Stops

ShapeOne® TAPS

Material Surgical steel

90.00 €			
	Ø 3.7 mm TAPS1B37	Ø 4.1 mm TAPS1B41	Ø 4.7 mm TAPS1B47

iMAX® TAPS

Material Surgical steel

90.00 €					
	Ø 3.3 mm TAPIMAX33	Ø 3.7 mm TAPIMAX37	Ø 4.1 mm TAPIMAX41	Ø 4.7 mm TAPIMAX47	Ø 5.2 mm TAPIMAX52

SHORT STOPS FOR DRILLS

Material Ti-Gr4

35.00 €						
	h 4 mm DS43341M	h 6 mm DS63341M	h 8 mm DS83341M	h 10 mm DS103341M	h 11 mm DS113341M	h 13 mm DS133341M





PARALLEL PINS

Material Ti-Gr4

20.00 €			
	0° PP	18° PP18	30° PP30

IMPLANT DRIVERS

Material Surgical steel

Internal Hex 2.1 mm Narrow  45.00 €		
	h 25 mm HDH21S	h 30 mm HDH21L
External Hex 2.4 mm Narrow  45.00 €		
	h 25 mm HDH24S	h 30 mm HDH24L

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps

Digital prosthetic  
solutions


Drills, countersinks  
taps, stops

Surgical  
kit

Surgical  
protocol

Package  
leaflets



Implant systems	Surgical Kit		Drills, Countersinks, Taps, Stops	
	description   €			
Internal hex prosthetic solutions	IMPLANT DRIVERS <i>(continued)</i>		Material Surgical steel	
External hex prosthetic solutions	For internal Hex 2.5 mm Regular 45.00 €	 h 25 mm HDH25S	 h 30 mm HDH25L	
Conical connection prosthetic solutions	For external Hex 2.7 mm Regular 45.00 €	 h 25 mm HDH27S	 h 30 mm HDH27L	
Internal octagon prosthetic solutions				
Shapemini prosthetic solutions	For octagon 3.1 mm 45.00 €	 h 25 mm HDH31S	 h 30 mm HDH31L	
MUA and OnePiece components				
Accessories, rings and retaining caps	For straight MUA 45.00 €	 HDH20		
Digital prosthetic solutions	For Mini implants 2 mm 45.00 €	 HDH25M		
Drills, countersinks taps, stops	For OnePiece da 4.5mm 45.00 €	 NHSMHDH		
Surgical kit				
Surgical protocol				
Package leaflets				

Surgical Kit

Drills, Countersinks, Taps, Stops

description | €

ORIENTER POSITION

Material Surgical steel

For OnePiece NHSM

21.00 €



NHSMFL

GUIDE TO DRILL INCLINATION

Material Surgical steel

0° - 18° - 30°

120.00 €




NHSMG

PROSTHETIC SCREWDRIVERS


Material Surgical steel

For Hex 1.25 mm contra-angle connection

45.00 €



short THDDS



long THDDL

For hex da 1.25 mm manual

45.00 €  
65.00 € (angled)



h 23 mm MDS



h 29 mm MDL



angled h 29 mm MDLA




30° angled MDLAA

SCREWDRIVERS

Material Surgical steel

Long


65.00 €



IDL

For iRETOR®

95.00 €



8393

Implant systems

Internal hex prosthetic solutions

External hex prosthetic solutions

Conical connection prosthetic solutions

Internal octagon prosthetic solutions

Shapemini prosthetic solutions

MUA and OnePiece components

Accessories, rings and retaining caps











Digital prosthetic solutions

Drills, countersinks taps, stops

Surgical kit

Surgical protocol

Package leaflets

Implant systems	Surgical Kit	Drills, Countersinks, Taps, Stops		
		description   €		
Internal hex prosthetic solutions	External hex prosthetic solutions	REMOVAL TOOLS		
		Material Surgical steel		
Conical connection prosthetic solutions	Internal octagon prosthetic solutions	95.00 €	 <p>for implants TRI</p>	 <p>For abutments TRT</p>
Shapemini prosthetic solutions	MUA and OnePiece components	EXTRACTION KIT FOR RETAINING SCREWS		
		Material Surgical steel		
Accessories, rings and retaining caps	Digital prosthetic solutions	Full kit 270.00 €	 <p>Drill D15RS 110.00 €</p>	 <p>Guide for drills GRS 60.00 €</p>
				 <p>Holder for GRS SGRS 100.00 €</p>
Drills, countersinks taps, stops	Surgical Kit	RATCHET		
		Material Surgical steel		
Surgical protocol	Package leaflets	Reverse fixed ratchet 90.00 €	 <p>R2</p>	
			 <p>WH2</p>	
Surgical protocol	Package leaflets	ratchet wrench 140.00 €	 <p>TW2</p>	
			 <p>RA2</p>	
Surgical protocol	Package leaflets	Complete ratchet 280.00 €	 <p>TWA2</p>	

PROSTHETIC KIT COMPOSITION



Code	Description
MDS	short screwdriver hex 1.25 mm h 23 mm
MDL	long screwdriver hex 1.25 mm h 29 mm
TRT	removal tool for abutments
HDH20	Implant driver for straight MUA
THDDS	short prosthetic screwdriver hex 1.25
THDDL	long prosthetic screwdriver hex 1.25
TWA2	ratchet wrench
MDLA*	angled screwdriver hex 1.25 mm h 29 mm

\*available on request

Implant systems

Internal hex prosthetic solutions

External hex prosthetic solutions

Conical connection prosthetic solutions

Internal octagon prosthetic solutions

Shapemini prosthetic solutions

MUA and OnePiece components

Accessories, rings and retaining caps

Digital prosthetic solutions

Drills, countersinks taps, stops

Surgical kit

Surgical protocol

Package leaflets

## ERGO KIT COMPOSITION



Code	Description
DE	drill extender
LD	lance drill
D20M	pilot drill Ø 2.0
D2428M	drill Ø 2.4 2.8 mm
D2833M	drill Ø 2.8 3.3 mm
D3338M	drill Ø 3.3 3.8 mm
D3844M	drill Ø 3.8 4.4 mm
D4448M	drill Ø 4.4 4.8 mm
CSD37	countersink Ø 3.7
CSD41	countersink Ø 4.1
CSD47	countersink Ø 4.7
CSD52	countersink Ø 5.2
TAPS1B37	tap Ø 3.7
TAPS1B41	tap Ø 4.1
TAPS1B47	tap Ø 4.7
TAPiMAX37	tap Ø 3.7
TAPiMAX41	tap Ø 4.1
TAPiMAX47	tap Ø 4.7
TAPiMAX52	tap Ø 5.2
DS43341M	stop 4.5
DS63341M	stop 6.5
DS83341M	stop 8
DS103341M	stop 10
DS113341M	stop 11.5
DS133341M	stop 13
TWA2	complete ratchet
THDDS	short screwdriver
THDDL	long screwdriver
HDH25S	short driver hex 2.5
HDH25L	long driver hex 2.5



SINUS LIFT COMPOSITION



Code	Description
DE	drill extender
SPD	standard pilot drill
SPI	standard start drill
SD5	standard drill h 5 mm
SD6	standard drill h 6 mm
SD7	standard drill h 7 mm
SD8	standard drill h 8 mm
SBL	standard body lift
APD	advanced pilot drill
AID	advanced start drill
AD2	advanced drill h 2 mm
AD3	advanced drill h 3 mm
AD4	advanced drill h 4 mm
ABL	advanced body lift
RBL	ratchet body lift

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps

Digital prosthetic  
solutions

Drills, countersinks  
taps, stops

Surgical  
kit

Surgical  
protocol

Package  
leaflets

heights from 8 to 16 mm

# SURGICAL PROTOCOL

ShapeOne

S1B - S1T - S1Tn

Ø	Heights
3.7	h 8 - 10 - 11.5 - 13 - 16 mm
4.1	
4.7	

iMAX

NHSl - NHSE

Ø	Heights
3.3	h 10 - 11.5 - 13 - 16 mm
3.7	h 8 - 10 - 11.5 - 13 - 16 mm
4.1	
4.7	
5.2	h 8 - 10 - 11.5 - 13 mm

iMAX

NHSlC

Ø	Heights
3.3	h 10 - 11.5 - 13 - 16 mm
3.7	h 8 - 10 - 11.5 - 13 - 16 mm
4.1	
4.7	
5.2	h 8 - 10 - 11.5 - 13 mm

iMAX MUA OnePiece

NHSM

Angle	Ø	Heights
0°	3.3	h 10 - 11.5 - 13 mm
	3.7	
	4.1	
18°	3.7	h 11.5 - 13 mm
	4.1	
30°	3.7	h 11.5 - 13 - 16 mm
	4.1	

Volution

SVB

Ø	Heights
3.3	h 10 - 11.5 - 13 - 16 mm
3.7	h 8 - 10 - 11.5 - 13 - 16 mm
4.1	
4.7	
5.2	h 8 - 10 - 11.5 - 13 mm

**DRIVERS- FINAL DRILLS- COUNTERSINKS - TAPS**

heights from 8 to 16 mm

Implant system	Drivers	Ø	Final drills (Bone d4)	Final drills (Bone d3-d2-d1)	CSD (Bone d3 - d2)	TAP (d1 bone)
iMAX NHSI 3.3	HDH21S HDH21L	3.3	D2024M	D2428M	CSD33	TAPIMAX33
SHAPEONE B	HDH25S HDH25L	3.7	D2428M	D2833M	CSD37	TAPS1B37
		4.1	D2833M	D3338M	CSD41	TAPS1B41
		4.7	D3338M	D3844M	CSD47	TAPS1B47
SHAPEONE Tn	HDH25S HDH25L	3.7	D2428M	D2833M	CSD37	TAPS1B37
		4.1	D2833M	D3338M	CSD41	TAPS1B41
		4.7	D3338M	D3844M	CSD47	TAPS1B47
iMAX NHSI	HDH25S HDH25L	3.7	D2428M	D2833M	CSD37	TAPIMAX37
		4.1	D2833M	D3338M	CSD41	TAPiMAX41
		4.7	D3338M	D3844M	CSD47	TAPiMAX47
		5.2	D3844M	D4448M	CSD52	TAPiMAX52
iMAX NHSIC Narrow	HDH21S HDH21L	3.3	D2024M	D2428M	CSD33	TAPIMAX33
		3.7	D2428M	D2833M	CSD37	TAPiMAX37
iMAX NHSIC Regular	HDH25S HDH25L	4.1	D2833M	D3338M	CSD41	TAPiMAX41
		4.7	D3338M	D3844M	CSD47	TAPiMAX47
		5.2	D3844M	D4448M	CSD52	TAPiMAX52
Volution SVB	HDH21S HDH21L	3.3	D2024M	D2428M	CSD33	
	HDH25S HDH25L	3.7	D2428M	D2833M	CSD37	
		4.1	D2833M	D3338M	CSD41	
		4.7	D3338M	D3844M	CSD47	
		5.2	D3844M	D4448M	CSD52	
iMAX NHSE 3.3	HDH24S HDH24L	3.3	D2024M	D2428M	CSD33	TAPIMAX33
iMAX NHSE	HDH27S HDH27L	3.7	D2428M	D2833M	CSD37	TAPiMAX37
		4.1	D2833M	D3338M	CSD41	TAPiMAX41
		4.7	D3338M	D3844M	CSD47	TAPiMAX47
		5.2	D3844M	D4448M	CSD52	TAPiMAX52
SHAPEONE T (abutment included in the pack)	HDH25S HDH25L	3.7	D2428M	D2833M	CSD37	TAPS1B37
		4.1	D2833M	D3338M	CSD41	TAPS1B41
		4.7	D3338M	D3844M	CSD47	TAPS1B47
SHAPEONE T (after removing the abutment tighten the implant)	HDH31S HDH31L	3.7	D2428M	D2833M	CSD37	TAPS1B37
		4.1	D2833M	D3338M	CSD41	TAPS1B41
		4.7	D3338M	D3844M	CSD47	TAPS1B47
iMAXMUA 0°	NHSMHDH	3.3	D2024M	D2428M	CSD33	TAPIMAX33
iMAXMUA 18°	NHSMFL	3.7	D2428M	D2833M	CSD37	TAPiMAX37
iMAXMUA 30°	(driver)	4.1	D2833M	D3338M	CSD41	TAPiMAX41
SHAPEMINI	HDH25M	2.7	D20M	D2024M		

 Implant  
systems

 Internal hex  
prosthetic solutions

 External hex  
prosthetic solutions

 Conical connection  
prosthetic solutions

 Internal octagon  
prosthetic solutions

 Shapemini  
prosthetic solutions

 MUA and OnePiece  
components

 Accessories, rings  
and retaining caps

 Digital prosthetic  
solutions

 Drills, countersinks  
taps, stops

 Surgical  
kit

 Surgical  
protocol

 Package  
leaflets


## SHAPEONE

Ø 3.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 3.7 countersink and 3.7 tap

For S1T implant system use the removal tool for abutments TRT after the implant has been placed in the bone



	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Countersink	Tap
Soft d4	•	•	•					
Medium d3,d2	•	•	•	•			•	
Compact d1	•	•	•	•			•	•


## SHAPEONE

Ø 4.1

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 4.1 countersink and 4.1 tap

For S1T implant system use the removal tool for abutments TRT after the implant has been placed in the bone



	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Countersink	Tap
Soft d4	•	•	•	•				
Medium d3,d2	•	•	•	•	•		•	
Compact d1	•	•	•	•	•		•	•


## SHAPEONE

Ø 4.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 4.7 countersink and 4.7 tap

For S1T implant system use the removal tool for abutments TRT after the implant has been placed in the bone








	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Countersink	Tap
Soft d4	•	•	•	•	•			
Medium d3,d2	•	•	•	•	•	•	•	
Compact d1	•	•	•	•	•	•	•	•

iMAX - iMAXMUA - iMAX NHSIC

Ø 3.3







Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 3.3 countersink and 3.3 tap

									
	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•							
Medium d3,d2	•	•	•					•	
Compact d1	•	•	•					•	•

iMAX - iMAXMUA - iMAX NHSIC

Ø 3.7








Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 3.7 countersink and 3.7 tap

									
	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•						
Medium d3,d2	•	•	•	•				•	
Compact d1	•	•	•	•				•	•

iMAX - iMAXMUA - iMAX NHSIC

Ø 4.1

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 4.1 countersink and 4.1 tap


									
	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•	•					
Medium d3,d2	•	•	•	•	•			•	
Compact d1	•	•	•	•	•			•	•



## iMAX - iMAX NHSIC

Ø 4.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 4.7 countersink and 4.7 tap




	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•	•	•				
Medium d3,d2	•	•	•	•	•	•		•	
Compact d1	•	•	•	•	•	•		•	•

## iMAX - iMAX NHSIC

Ø 5.2

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 5.2 countersink and 5.2 tap







	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•	•	•	•			
Medium d3,d2	•	•	•	•	•	•	•	•	
Compact d1	•	•	•	•	•	•	•	•	•

VOLUTION

Ø 3.3






Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 3.3 countersink

								
	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•						
Medium d3,d2	•	•	•					•
Compact d1	•	•	•					•

VOLUTION

Ø 3.7







Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 3.7 countersink

								
	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•					
Medium d3,d2	•	•	•	•				•
Compact d1	•	•	•	•				•

VOLUTION

Ø 4.1


Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 4.1 countersink

								
	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•	•				
Medium d3,d2	•	•	•	•	•			•
Compact d1	•	•	•	•	•			•

## VOLUTION

Ø 4.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 4.7 countersink




	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•	•	•			
Medium d3,d2	•	•	•	•	•	•		•
Compact d1	•	•	•	•	•	•		•

## VOLUTION

Ø 5.2

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 5.2 countersink



	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•	•	•	•		
Medium d3,d2	•	•	•	•	•	•	•	•
Compact d1	•	•	•	•	•	•	•	•

ShapeOne		S1B
Ø	Heights	
4.1	h 6.5 mm	
4.7		

ShapeOne		S1T - S1Tn
Ø	Heights	
4.1	h 4.5 - 6.5 mm	
4.7		

iMAX		NHSI - NHSE
Ø	Heights	
3.7	h 6.5 mm (NHSE only)	
4.1	h 6.5 mm	
4.7		
5.2		

iMAX		NHSIC
Ø	Heights	
4.1	h 6.5 mm	
4.7		
5.2		

Volution		SVB
Ø	Heights	
3.3	h 6.5 mm	
3.7		
4.1		
4.7		
5.2		

## DRIVERS- FINAL DRILLS- COUNTERSINKS - TAPS

short implants

Implant system	Drivers	Ø	Final drills (Bone d4)	Final drills (Bone d3-d2-d1)	CSD (Bone d3 - d2)	TAP (Bone d1)
SHAPEONE B	HDH25S	4.1	D3338M	D3844M	CSD41	TAPS1B41
	HDH25L	4.7	D3844M	D4448M	CSD47	TAPS1B47
SHAPEONE Tn	HDH25S	4.1	D3338M	D3844M	CSD41	TAPS1B41
	HDH25L	4.7	D3844M	D4448M	CSD47	TAPS1B47
iMAX NHSI	HDH25S HDH25L	3.7	D2833M	D3338M	CSD37	TAPiMAX37
		4.1	D3338M	D3844M	CSD41	TAPiMAX41
		4.7	D3844M	D4448M	CSD47	TAPiMAX47
		5.2	D4448M	D4448M	CSD52	TAPiMAX52
iMAX NHSIC Regular	HDH25S HDH25L	4.1	D3338M	D3844M	CSD41	TAPiMAX41
		4.7	D3844M	D4448M	CSD47	TAPiMAX47
		5.2	D4448M	D4448M	CSD52	TAPiMAX52
Volution SVB	HDH21S HDH21L	3.3	D2428M	D2833M	CSD33	
		3.7	D2833M	D3338M	CSD37	
	HDH25S HDH25L	4.1	D3338M	D3844M	CSD41	
		4.7	D3844M	D4448M	CSD47	
iMAX NHSE	HDH27S HDH27L	5.2	D4448M	D4448M	CSD52	
		3.7	D2833M	D3338M	CSD37	TAPiMAX37
		4.1	D3338M	D3844M	CSD41	TAPiMAX41
		4.7	D3844M	D4448M	CSD47	TAPiMAX47
SHAPEONE T (abutment included in the pack)	HDH25S HDH25L	5.2	D4448M	D4448M	CSD52	TAPiMAX52
		4.1	D3338M	D3844M	CSD41	TAPS1B41
SHAPEONE T (after removing the abutment tighten the implant)	HDH31S HDH31L	4.7	D3844M	D4448M	CSD47	TAPS1B47
		4.1	D3338M	D3844M	CSD41	TAPS1B41



# SHAPEONE










short implants

Ø 4.1

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 4.1 countersink and 4.1 tap

For S1T implant system use the removal tool for abutments TRT after the implant has been placed in the bone

													
Bone	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9				
Soft d4	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap				
Soft d4	•	•	•		•								
Medium d3,d2	•	•	•		•	•		•					
Compact d1	•	•	•		•	•		•		•			•

# SHAPEONE










short implants

Ø 4.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 4.7 countersink and 4.7 tap

For S1T implant system use the removal tool for abutments TRT after the implant has been placed in the bone

											
Bone	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9		
Soft d4	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap		
Soft d4	•	•	•	•		•					
Medium d3,d2	•	•	•	•		•	•	•			
Compact d1	•	•	•	•		•	•	•	•		•










# iMAX

short implants

Ø 3.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 3.7 countersink and 3.7 tap

													
Bone	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9				
Soft d4	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap				
Soft d4	•	•		•									
Medium d3,d2	•	•		•	•			•					
Compact d1	•	•		•	•			•		•			•


## iMAX - iMAX NHSIC

short implants

Ø 4.1

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 4.1 countersink and 4.1 tap



	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•		•				
Medium d3,d2	•	•	•		•	•		•	
Compact d1	•	•	•		•	•		•	•


## iMAX - iMAX NHSIC

short implants

Ø 4.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 4.7 countersink and 4.7 tap



	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•	•		•			
Medium d3,d2	•	•	•	•		•	•	•	
Compact d1	•	•	•	•		•	•	•	•


## iMAX - iMAX NHSIC

short implants

Ø 5.2

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone

Use 5.2 countersink and 5.2 tap











	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink	Tap
Soft d4	•	•	•	•	•		•		
Medium d3,d2	•	•	•	•	•		•	•	
Compact d1	•	•	•	•	•		•	•	•

VOLUTION

short implants

Ø 3.3

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 3.3 countersink









								
Bone	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•		•					
Medium d3,d2	•		•	•				•
Compact d1	•		•	•				•

VOLUTION

short implants

Ø 3.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 3.7 countersink









								
Bone	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•		•				
Medium d3,d2	•	•		•	•			•
Compact d1	•	•		•	•			•

VOLUTION

short implants

Ø 4.1

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 4.1 countersink

								
Bone	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•		•			
Medium d3,d2	•	•	•		•	•		•
Compact d1	•	•	•		•	•		•

VOLUTION

short implants

Ø 4.7

Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 4.7 countersink

	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•	•		•		
Medium d3,d2	•	•	•	•		•	•	•
Compact d1	•	•	•	•		•	•	•

VOLUTION

short implants

Ø 5.2


Sink countersink: up to 1.4mm for d3 medium bone/ up to 2.8mm for d2 medium bone and d1 compact bone  
Use 5.2 countersink

	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8
Bone	Lance drill	Drill 2.0 2.4	Drill 2.4 2.8	Drill 2.8 3.3	Drill 3.3 3.8	Drill 3.8 4.4	Drill 4.4 4.8	Countersink
Soft d4	•	•	•	•		•		
Medium d3,d2	•	•	•	•		•	•	•
Compact d1	•	•	•	•		•	•	•

MINIMUM IMPLANTS SIZE ALLOWED FOR POSITION

iRES declines all responsibility in case of failure if the information leaflet on page 73, are not be respected as regard the implants position in relation to implants site and diameters

..... Implant head



Implant systems
Internal hex prosthetic solutions
External hex prosthetic solutions
Conical connection prosthetic solutions
Internal octagon prosthetic solutions
Shapemini prosthetic solutions
MUA and OnePiece components
Accessories, rings and retaining caps
Digital prosthetic solutions
Drills, countersinks taps, stops
Surgical kit
Surgical protocol
Package leaflets



A black and white photograph of a man with glasses and a beard, looking down at a smartphone he is holding in his right hand. He is wearing a checkered shirt. The background is blurred, showing what appears to be a modern building with large windows.

## NEW iRES® eShop

24/7

---

[www.eshop.ires.dental](http://www.eshop.ires.dental)

Sign up now and find out  
**all the benefits of online**  
shopping. You can easily  
order all our products thanks  
to a new clear, simple and  
functional graphical interface

---





The SHAPE1, IMAX, VOLUTION, iMAXMUA and SHAPEM- INI implant systems, by I-RES Sagl, offers the dentist a wide choice of titanium implant configurations that differ in diameter, height and possibility of surgical po- sitioning A) submerged/bone level, B) transmucosal/ tissue level, and various prosthetic components for the different rehabilitation procedures

The sterility of the implant is guaranteed by the se-

Implant surgery operations may lead to some complications that are usually temporary and restricted to the area of operation, such as inflammation, paresthesia, haematoma; there may also be injuries to nerves, to vascular complexes and the membrane of the maxillary sinus. Bone sequestration has rarely

Ires' Sagl implants are DISPOSABLE devices. Their re-use is not desirable from a medical, legal and ethical point of view. The use of not validated sterilization procedures can cause both the infection onset in the patient and impair the product performances. The failure compliance with these instructions implies a different use as provided by the manufacturer and those who make the reuse must take this action on their own responsibility.



Do not use if the packaging is damaged

The I-RES implant-prosthetic system is intended for use in the oral cavity and is provided with useful components to enable the dentist and the dental technician to prosthetically complete the operation begun by the dentist by inserting the implant in the patient. The I-RES implant-prosthetic system completes the line of I-RES dental implants and the respective I-RES instruments.

#### Product characteristics

Healing screws. The healing screw is a device used by the dentist to keep the oral mucosa pervious in the vicinity of the implant previously inserted. Once the soft tissues have healed, this will allow the dentist to perform the subsequent manoeuvres for prosthetization. The healing screws are made of grade 5 titanium.

Transfer. The transfer in grade 5 titanium is the instrument that allows the transfer from the mouth to a model of the information needed for the prosthetic connection and for making the respective prosthesis. There are two types of transfer: "closed tray and open tray", and they are all composed of two parts (a screw and a repositioner). After being inserted in the implant and secured to it with the screw, the transfer is ready to take the impression in the mouth.

Analog. The analog made of grade 5 titanium has the function of reproducing the internal characteristics of the implant and it must be securely fixed to the transfer. Once joined, the model can be cast.

Straight, angled and millable posts. They are made of grade 5 titanium; they have different shapes depending on the characteristics they have to satisfy, they are used mostly for prosthetic rehabilitations of bridges or crowns. The choice of the device that must be connected to the analog in the first phase is dictated by the clinical and processing decisions, which are at the discretion of the dentist and the dental technician.

Plastic posts. Plastic posts may be divided into two families, one for using directly in the oral cavity, appropriately modified and connected to the post to support temporary prostheses, one for the transformation of plastic posts into metal posts by the dental technician, with processing characteristics that are at the discretion of the dentist and the dental technician.

Gold Bases. These are components made of gold alloy and allow the creation of customized posts using overcasting techniques.

Ball attachments. Ball attachments are made of grade 5 titanium and, once fixed to the implants, they are able to act as an anchorage by means of special attachments to the patient's mobile prosthesis.

#### Contraindications:

Do not use I-RES products on patients who have allergies to the materials of which the component is made. The use of I-RES components in patients who have metabolic and periodontal diseases or poor oral hygiene may be prejudicial to the success of the

work, as may prosthetic constructions not in line with international standards. The lack of periodic controls, which the patient must undergo with his dentist after prosthetisation, may compromise the life of the implant-prosthetic system.

#### Warnings:

I-RES prosthetic components are reserved for use by personnel with knowledge of the subject. I-RES points out that alterations to the implant/post connections may be prejudicial to the success of the work, as may the failure to use original components. When using prosthetic components it is important to follow the instructions given by the dentist and the dental technician. When using prosthetic components

in the oral cavity it is important to respect the final tightening value which must be between 20 and 30 Ncm, as better specified in the catalogue.

#### Collateral effects

Today there are no known collateral effects in the use of I-RES components that can endanger the patient's health.

#### Prosthetic planning:

The choice of the I-RES components to be used for the case is the specific responsibility of the dentist and of the dental technician, depending on their requirements.

#### Materials and packaging:

All I-RES prosthetic components are packed in such a way as to be immediately identifiable, once removed from their pack; it is important for the operator to pay great attention in identifying them to avoid changes of position during work. It is useful to make note of the material batch used on the patient's file, for the purpose of traceability. Whether it has been processed or not, before inserting the I-RES prosthetic component in the oral cavity it is of fundamental importance that it be washed and sterilized. Some I-RES components are single-use, so intended for only one patient.

#### Cleaning | sterilization | storage:

**Caution !!! All prosthetic components for dental implants are sold NON-STERILE.**

Before use, all prosthetic components must be cleaned, disinfected and sterilized. These processes must also be performed before intraoral use, i.e. before each use for any test phases and in any case before final restoration loading. Repetition of the processes described in this paragraph does not alter the characteristics of these devices. Failure to follow these indications may lead to the onset of infections and complications for the implant and, more generally, for the patient. Important: care must be taken during the subsequent phases in preserving the zone of the connection with the implant (hexagon/octagon/threading).

#### a. Cleaning:

In case of automatic cleaning: use an ultrasound bath with a suitable detergent solution. Use neutral detergents only. Follow the manufacturer's instructions concerning concentrations and washing times. Use demineralised water to prevent the formation of stains and marks.

When cleaning manually: use a suitable neutral detergent and follow the manufacturer's user instructions. Brush the products with a soft-bristled brush (non-metal bristles) under running water. Use the brush to apply the detergent to all surfaces. Rinse with distilled water. After rinsing, dry the devices thoroughly and place them inside suitable sterilization bags.

#### b. Sterilization:

Place in a vacuum autoclave and sterilize as follows: Temperature = 121 – 124°C, with autoclave cycle of at least 20 minutes and drying cycle of 15 minutes.

#### c. Storage:

After sterilization, the product must remain in the sterilization bags. Only open the bags immediately prior to use. In normal conditions, sterilization bags maintain the sterility of the contents, unless the wrapping is damaged. Therefore, do not use components if the bags in which they were kept are damaged, and re-sterilizes in new bags before using them again. The storage time of products sterilized inside the bags should not exceed that recommended by the manufacturer of the bags.

The product must be stored in a cool dry place, away from direct sunlight, water and heat sources.

#### ATTENTION:

Some components such as transfers and healing screws are devices that can be reused after.


#### CLEANING/STERILIZATION/STORAGE (follow the re-pective indications).

**DO NOT REUSE a device classified as SINGLE-USE. Although it cannot be seen, it could be mechanically deformed or have been contaminated.**

#### Disposal procedures:

If removed from the oral cavity due to biological or mechanical failure, the prosthetic components must be disposed of as biological waste according to local regulations. More detailed information on the use of the medical device can be found in the specific Surgical Protocol available on the site [www.i-res-group.com](http://www.i-res-group.com) or in the IRES Shape1 catalogue supplied by the Manufacturer.

#### Symbols on the package:

	MANUFACTURER I-RES® SAGL Riva Caccia, 1/D 6900 Lugano [Switzerland] <a href="mailto:info@ires.dental">info@ires.dental</a> <a href="http://www.ires.dental">www.ires.dental</a>
	CE Mark according to standard MDD93/42/EEC
	Batch number
	use before the expiry date
	Do not reuse
	Follow the instructions given in the illustrative leaflet
	Do not expose to direct sunlight
	Do not expose to rain and keep in an environment free from damp
	Do not use if the packaging is damaged
	not sterile

**INSTRUCTIONS FOR IRES ROTARY INSTRUMENTS (DRILLS - COUNTERSINKS - TAPS) FOR THE PREPARATION OF THE SITE THAT HAS TO RECEIVE IRES® SHAPE1® IMPLANTS**

**Product description:**

Dental drills, produced by I-RES Sagl, must be used as tools to perforate the bone. The diameters to be used, the lengths and the drilling sequence (number of drills to be used) are the sole choice and decision of the dentist, depending on the surgical protocol that must be followed. The maximum recommended speed is 800 rpm with saline solution applied directly on the drill to assist cooling.

a) The sole purpose of the initial precision drill is to incise the cortical bone in a very precise point where it will later be drilled.

b) The helical drills have laser markings for reference which identify the depth to be reached. Of course, in the use of this type of drill the manual skill and experience of the dental surgeon are extremely important, especially for stopping at the chosen depth.

c) Countersinks are used when it is necessary to widen the initial part of the hole made to adapt the shape that of the neck of the implant to be inserted. The maximum recommended speed is 300 rpm with saline solution applied directly on the drill to assist cooling.

d) Bone taps: in particularly dense bone (type I), before insertion it is advisable to use a bone tap with the same profile as the implant to be inserted. The bone tap has a greater cutting power than the implant, allowing the site to be prepared with reduced trauma. The maximum recommended speed is 30 rpm with saline solution applied directly on the bone tap to assist cooling.

**Materials used:**

All I-RES Sagl drills are made of medical grade steel and undergo hardening heat treatment. The maximum recommended number of uses of the devices is 40 times.

**Warnings and general precautions:**

- It is fundamental to respect the surgical protocol that establishes the diameters, lengths and the sequence of use. The operator is fully responsible for any uses other than those indicated.

- Check that the drills to be used are in good condition, already cleaned and sterilized.

- Check that the drills are in good condition and have not been used more than 40 times.

- Before using them, check that the hand-piece holds the drills perfectly secure and that they rotate in the correct direction.

- Ensure that there is adequate irrigation.

- The application of leverage during drilling could cause breakage of the drill, the hand-piece, or the

bone walls on which you are working.

During drilling always exert alternating pressure, using the intermittent drilling technique.

- Always check that the laser marking that indicates diameter and length is clearly visible.

- Any eccentricity or lack of straightness in the drill could result in an oversized hole.

- Wear eye protection, to protect against particles that may be ejected.

**CLEANING / STERILIZATION / STORAGE:**

The medical devices are supplied NON-STERILE.

Before use, all rotary devices must be cleaned, disinfected and sterilized.

Failure to follow these indications may lead to the onset of infections and complications for the implant and, more generally, for the patient.

**a. Cleaning**

In case of automatic cleaning: use an ultrasound bath with a suitable detergent solution. Use neutral detergents only. Follow the manufacturer's instructions concerning concentrations and washing times. Use demineralised water to prevent the formation of stains and marks.

When cleaning manually: use a suitable neutral detergent and follow the manufacturer's user instructions. Brush the products with a soft-bristled brush (non-metal bristles) under running water. Use the brush to apply the detergent to all surfaces. Rinse with distilled water. After rinsing, dry the devices thoroughly and place them inside suitable sterilization bags.

**b. Sterilization**








Place in a vacuum autoclave and sterilize as follows: Temperature = 121 – 124°C, with autoclave cycle of at least 20 minutes and drying cycle of 15 minutes.

**c. Storage**

After sterilization, the product must remain in the sterilization bags. Only open the bags immediately prior to use. In normal conditions, sterilization bags maintain the sterility of the contents, unless the wrapping is damaged. Therefore, do not use components if the bags in which they were kept are damaged, and re-sterilize in new bags before using them again. The storage time of products sterilized inside the bags should not exceed that recommended by the manufacturer of the bags. The product must be stored in a cool dry place, away from direct sunlight, water and heat sources.

More detailed information on the use of the medical device can be found in the Surgical Protocol. If you do not have a copy, request one from your distributor or directly from the manufacturer.

**Symbols on the package:**

	MANUFACTURER I-RES® SAGL Riva Caccia, 1/D 6900 Lugano [Switzerland] info@ires.dental www.ires.dental
	CE Mark according to standard MDD93/42/EEC
	Batch number
	Follow the instructions given in the illustrative leaflet
	Do not expose to direct sunlight
	Do not expose to rain and keep in an environment free from damp
	Do not use if the packaging is damaged
	Not sterile

Implant  
systems

Internal hex  
prosthetic solutions

External hex  
prosthetic solutions

Conical connection  
prosthetic solutions

Internal octagon  
prosthetic solutions

Shapemini  
prosthetic solutions

MUA and OnePiece  
components

Accessories, rings  
and retaining caps

Digital prosthetic  
solutions

Drills, countersinks  
taps, stops

Surgical  
kit

Surgical  
protocol

Package  
leaflets





## CERTIFICATION

iRES® SAGL is a company with system EN ISO 13485 quality management whose quality standards are rigorously tested by the Notified Body. Process of on going FDA certification.











implant  
systems



prosthetic  
solutions



continuing  
education



guided  
surgery



regenerative  
materials



laser  
fixlite®

**info@ires.dental**  
**www.ires.dental**  
**eshop.ires.dental**

**ires**  | complete  
solutions  
for oral  
surgery

