

成型尺寸: 100×155mm 材质:80g双胶纸 黑白印刷 骑钉成本

物料名称	BCL-A英文说明书	设计	日期	20190924
物料编号	PA-BCL.A-M-E-1	审核	日期	
版本	A1	批准	日期	

Contents

1 Symbol instruction 1
2 Introduction 2
3 Intended use, structure and contraindication2
4 Technical specifications 3
5 Components 3
6 Install and uninstall way 4
7 Cleaning and sterilization 4
8 Operation 5
9 Precaution 6
10 Maintenance 7
11 Trouble shooting 7
12 Storage and transportation 7
13 After-sale service 8
14 Environment protection 8
15 Manufacturer' s right8
16 EMC-Declaration of conformity8

Ι

b Trade	mark		Caution
Type B	applied		Follow instructions for use
Use in	door only	IPX0	Ordinary equipment
135°C	avable	SP.	Screw inside/outside
Date of r	manufacture		Manufacturer
*			y be harmful to the eyes ght source
Ŕ	Appliance WEEE dir	complia ective	nce with
70106	Atmosphe 70kPa-10		sure for storage
80%	Humidity 0% [~] 80%	limitatio	n for storage
-10℃	Temperatu -10℃~+5		tion for storage



4. Technical specifications 4.1 Power input : AC24V 50Hz. 4.2 Input power: 10VA. 4.3 Dimensions : 26×221mm. 4.4 Net weight : 0.18kg. 4.5 Components: details in packing list. 4.6 LED light performance: 4.6.1 LED power: 5W. 4.6.2 Wave length : 440nm~480nm. 4.6.3 Category : I. 4.6.4 Light intensity: a) Low power:1000mW/cm²~1300mW/cm². b) High power:1600mW/cm²~1800mW/cm². 4.6.5 Optical active area 56mm². 4.7 Working condition: 4.7.1 Environment temperature : $+5^{\circ}C^{+30}$. 4.7.2 Relative humidity : 0~80%. 4.7.3 Atmosphere pressure : 70KPa~106KPa. 4.8 Design safety category: 4.8.1 Running mode : Indirect running 4.8.2 Protection type against electrical shock:class I. 4.8.3 Protection against electrical shock:type B. 4.8.4 Protection against harmful ingress of water or particular matter:ordinary equipment (IPX0). 4.8.5 Safety in the presence of flammable anesthetic mixture with air,oxygen or nitrous oxide:Non AP/APG type equipment. 5.Components 5.1 No spare parts, accessories, consumables. 5.2 More details about product and accompanying documents in the packing list. 3



8.Operation

8.1 Open the packing and check if all the parts and accessories are complete according to the packing list. 8.2 Ensure protection measures before use.For example, wear anti-glare glasses. 8.3 Pull out of the lamp cap and pack it in sterilization pouch or sterilization box to sterilize it in autoclave. 8.4 Install the lamp cap on the main unit. 8.5 Press 'M' button and choose high power mode or low power mode. 8.6 Press 'T' button and choose time setting. 8.7 Put the light source alignment resin which needs to be solid, press 'P' button, and then there is a 'Di' voice. The main unit shines blue light and starts countdown. 8.8 It can be stopped working at any time by short pressing 'P' button while operating. 8.9 The next working cyble should be started after a working

cycle ends for 15 seconds in order to avoid discomfort or scald to the patient due to the overheated lamp cap.

8.10 The lamp cap gets a little hot when working normally.But if the main unit or the lamp cap gets hot obviously, please turn off the device until the main unit becomes cool.8.11 If there is resin on the lamp cap after using,please

clean it by cotton cloth to make sure the power intensity is good.

8.12 After using curing light every time, pull out of the lamp cap, clean it and then pack it in sterilization pouch or sterilization box to sterilize it in case using it next time.

5

9.Safety precaution

GB9706.1 standards.

9.1 The lamp cap is an electrical product. Wipe it by cotton cloth or medical alcohol. Dipping in liquid is forbidden.
9.2 There is high nonopaque quartz glass at the front end of lamp cap. Knocking,striking or dropping it is forbidden.
9.3 During operation the light should be aimed straightly on the resin, to ensure solidification effectively.

9.4 Please wear anti-glare glasses to avoid blue light hurts eyes. Strictly prohibit aiming light at eyes directly.
9.5 Dedicated ~24V connector of dental unit supplies electricity to BCL-A built-in curing light. Power supply source of dental unit should conform to requirement of

9.6 Using product must conform to the related operation rules issued by health sector and the related laws requirements. It can be used by trained doctors or technicians only.

9.7 Suggest inspect light intensity of BCL-A curing light with light meter every month. If light intensity is not in the scope of setting, please send it back to us for repair.

9.8 Please use BCL-A curing light in the scope of cable length. Don't drag it please. Put it back to the hanger of dental unit after use.

9.9 Be sure to replace the original lamp cap made by us. If not original or not unmatched one, maybe it will damage or cause insufficient light intensity.

9.10 If necessary, circuit diagram, components and parts list, drawing statement, correction detailed rules and regulations can be offered to help technicians to repair product.



10.Maintenance

10.1 This product doesn't include repair parts. Machine maintenance and repair shall be made by the designated professional or repaired at a repair shop.

10.2 Users can replace lamp cap by themselves, but ensure to use original accessories.

10.3 Wipe surface of the lamp cap by purified water or disinfectant.

11.Trouble shooting

Fault	Possible cause	Solution
No	BCL-A curing light without electric power.	Check if BCL-A curing light connects dental unit correctly or not.
indication Non-act.	Curing light works for a long time and automatic protection starts.	Curing light stops working for several minutes,and it will return to normal.
	BCL-A curing light damages.	Contact the local dealer or us.
Light	LED ageing.	Replace lamp cap.
intensity insufficient	There is residual foreign matter at the front end of lamp cap.	Remove foreign matter.

If the problem still can't be solved, please contact the local dealer or us.

12 Storage and transportation

12.1 Environmental conditions of storage and transportation: a) Relative humidity:0% to 80%.

- b) Atmospheric pressure:70kPa to 106kPa.
- c) Environmental temperature:-10°C to +50°C.

12.2 Prevent excessive shock and vibration in transportation.

be sure to handle with care and avoid inversion.

12.3 Don't mix with dangerous goods during transportation.

12.4 Avoid the sun, rain or snow during transportation.

12.5 The equipment should be handled carefully and lightly.Be sure that it's far from the vibration, and install or store in a cool,dry and ventilated place.

12.6 Don't store the machine with the articles that are combustible, poisonous, caustic and explosive.

13.After-sale service

From the date of this equipment sold, if the equipment can't work normally because of the quality problems, our company will be responsible for the maintenance base on the warranty card. Please refer the warranty scope and period from warranty card.

14.Environmental protection

The equipment has no any harmful ingredient. Please dispose according to the local laws.

15.Manufacturer's right

We reserve the rights to change the design of the equipment, the technique, fittings, the instruction manual and the content of the original packing list at any time without notice. If there are some differences between blueprint and real equipment, take the real equipment as the norm.

16.EMC-Declaration of conformity Notes:

1) Maybe it will cause electromagnetic compatibility problem for this equipment or other equipments if change this equipment without our authorization.

2) Design and testing of this equipment conforms to the related operating instruction of electromagnetic compatibility.3) This equipment should not be used near other equipments or stacked up with other equipments. If needed to use the equipment under situations mentioned above, the user should verify that the equipment can work normally.

16. 1 Electric cable information

No	Description	Electric cable length (m)	Shielding (Yes/No)	Remarks
1	Power line (input end)	2. 0 m	No	None

8

Performance		Description	Response time of achieving normal performance
Light output	Press "P" light shine of lamp button a operation Set time choose hi	unit is turned on. button,then strong s from the front end cap. Pressing "P" gain can cancel when there is light. by 'T' button and igh intensity or low by 'M' button when b light.	<1s
ccessories,pc	apacitance.If	f use or change non-mat r and capacitance,that w	ill significantly dec
nductor and ca accessories,po electromagneti quipment par 6.4 Guidance Guidance an The model BCL	apacitance.If ower inductor ic emissions ts arbitrarily. and manufa d manufacture -A is intended	r and capacitance,that w and immunity performa	rill significantly dec nce.Please not cha stromagnetic emiss netic emissions letic environment
nductor and ca accessories,po electromagneti quipment par 6.4 Guidance Guidance an The model BCL specified below that it is used in	apacitance. If wer inductor c emissions ts arbitrarily. and manufa d manufacture -A is intended. The custome such an envi	r and capacitance, that w and immunity performan acturer's declaration-electormag d for use in the electromagr r or the user of the model E ronment.	rill significantly dec nce.Please not cha etromagnetic emiss netic emissions letic environment ICL-A should assure
nductor and ca accessories,po electromagneti quipment par 6.4 Guidance Guidance an The model BCL specified below	apacitance.If ower inducto c emissions ts arbitrarily. and manufa d manufacture -A is intended. The custome	r and capacitance, that w and immunity performan acturer's declaration-elector er's declaration-electromag f for use in the electromagr r or the user of the model E	rill significantly dec nce.Please not cha etromagnetic emissions letic environment BCL-A should assure energy only for its its RF emissions are o cause any
nductor and ca accessories,pc electromagneti quipment par 6.4 Guidance Guidance an The model BCL specified below that it is used in Emissions test RF emissions	apacitance. If wer inductor c emissions ts arbitrarily. and manufa d manufacture -A is intended. The custome such an envi Compliance	r and capacitance,that w and immunity performan acturer's declaration-elect er's declaration-electromag d for use in the electromagr r or the user of the model B ronment. Electromagnetic enviro The model BCL-A uses RF internal function. Therefore very low and are not likely j	rill significantly dec nce.Please not cha etromagnetic emissions letic emissions letic environment ICL-A should assure nment - guidance energy only for its its RF emissions are o cause any ronic equipment. It use in all estic and those directly tage power supplt
nductor and ca accessories,pc electromagneti quipment par 6.4 Guidance Guidance an The model BCL specified below that it is used in Emissions test RF emissions CISPR 11 RF emissions	apacitance. If ower inductor c emissions ts arbitrarily. and manufa d manufacture -A is intended .The custome such an envi Compliance Group 1	r and capacitance, that w and immunity performan acturer's declaration-elector er's declaration-electromage for use in the electromage r or the user of the model E ronment. Electromagnetic enviro The model BCL-A uses RF internal function. Therefore, very low and are not likely interference in nearby elect The model BCL-A is suitable for establishments other than dom connected to the public low-vo network that supplies buildinge	rill significantly dec nce.Please not cha etromagnetic emissions letic emissions letic environment ICL-A should assure nment - guidance energy only for its its RF emissions are o cause any ronic equipment. It use in all estic and those directly tage power supplt

	Guidance & Decla	aration — electr	omagnetic immunity
			agnetic environment specified beliw. assure that It is used in such an
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.		Not applicable	Not use mains power.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m ,3 A/m	3 A/m ,3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a	.c. mains voltage p	prior to application	n of the test level.
		10	

	Guidance & De	claration - Elec	tromagnetic immunity
			ctromagnetic environment specified below. hould assure that it is used in such an
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz 6 Vrms in ISM and amateur	3 Vrms 150 kHz to 80MHz 6 Vrms in ISM and amateur	Portable and mobile RF communications equipment should be used no closer to any part of the model BCL-A including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the
	radio bands	radio bands	transmitter.
	3 V/m	3 V/m	Recommended separation distance d=[3,5/V1]×P ^{1/2}
Radiated RF	80 MHZ to 2.5	80 MHZ to 2.5	d=1.2×P ^{1/2} 80 MHz to 800 MHz
IEC 61000-4-3	GHz 10 V/m	GHz 10 V/m	<i>d</i> =2.3×P ^{1/2} 800 MHz to 2.7 GHz
	80 MHZ to 2.7 GHz	80 MHZ to 2.7 GHz	where P is the maximum output power rating of the transmitter In watts (W) according to the transmitter manufacturer
	385MHz-	385MHz-	and d is the recommended separation distance in meters (m).
	5785MHz Test specifications for ENCLOSURE PORT	5785MHz Test specifications for ENCLOSURE PORT	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range.
	IMMUNITY to RF wireless communication	IMMUNITY to RF wireless communication	Interference may occur In the vicinity of equipment marked with the following symbol:
	equipment (Refer to table 9 of IEC 60601- 1-2:2014)	equipment (Refer to table 9 of IEC 60601- 1-2:2014)	(((•)))

Г

NOTE 2 Th	0 MHz end 800 MHz. the higher frequency range applies. se guidelines may not apply in all situations. Electromagnetic propagation is absorption and reflection from structures, objects and people.
telephone broadcast environma considere exceeds t to verify n be necess	ths from fixed transmitters, such as base stations for radio (cellular/cordless) and land mobile radios, amateur radio, AM and FM radio broadcast and TV cannot be predicted theoretically with accuracy. To assess the electromagnetic nt due to fixed RF transmitters, an electromagnetic site survey should be . If the measured field strength in the location in which the model BCL-A is used e applicable RF compliance level above, the models BCL-A should be observe etr mal operation. If abnormal performance is observed, additional measures may ary, such as reorienting or relocating the model BCL-A. equency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.
	12

Recomm portable and mobile RF	ended separation d		
The model BCL-A is inten- radiated RF disturbances BCL-A can help prevent e distance between portable	ded for use in electro is controlled. The cus lectromagnetic interfore and mobile RF com	magnetic env stomer or the erence by mai munications e	ironment in which user of the model ntaining a minimum quipment (transmitters)
and the model BCL-A as r power of the communication		according to	the maximum output
Rated maximum output power			ding to frequency
of transmitter W	150kHz to 80MHz d=1.2×P ^{1/2}	80MHz to 800MHz d=1.2×P ^{1/2}	800MHz to 2,5GHz d=2.3×P ^{1/2}
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10			
10	3.8	3.8	7.3
100 For transmitters rated at a is separation distance <i>d</i> in me the frequency of the transmitter transmitter in watts (W) acc NOTE I At 80 MHz and 800	12 maximum output pow eters (m) can be estii nitter, where <i>P</i> is the cordable to the trans 0 MHz. the separation	12 ver not listed a mated using th maximum out mitter manufa n distance for	23 bove, the recommended re equation applicable to put power rating of the cturer. the higher frequency
	12 maximum output pow eters (m) can be estin nitter, where <i>P</i> is the cordable to the trans 0 MHz. the separation asse guidelines may no in is affected by	12 wer not listed a mated using th maximum our mitter manufa n distance for ot apply in all	23 bove, the recommended he equation applicable to tput power rating of the cturer. the higher frequency situations.
100 For transmitters rated at a separation distance <i>d</i> in m the frequency of the transmitter transmitter in watts (W) acc NOTE I At 80 MHz and 800 range applies. NOTE 2 The Electromagnetic propagation	12 maximum output pow eters (m) can be estin nitter, where <i>P</i> is the cordable to the trans 0 MHz. the separation asse guidelines may no in is affected by	12 wer not listed a mated using th maximum our mitter manufa n distance for ot apply in all	23 bove, the recommended he equation applicable to tput power rating of the cturer. the higher frequency situations.

