

MAIN BUSINESS

OCEANREEF underwater, naturally.

The diving equipments produced by Mestel Safety are distributed under the Registered Trademark, OCEAN REEF.

5 A F E T Y electronics

- Electronics
- Under Water Communication Systems

SAFETY engineering

- Engineering
- R&D in Material Science and Tecnology
- Project Management

S A F E T Y protection equipment

The Military and Respiratory PPEs are known under the Mestel Safety brand or, simply, as SAFETY Protection Equipments.



EXPERT IN FULL FACE MASK







SAFETY FIRST, AS THE INNOVATOR OF THE MARKET

CO₂ LEVELS TESTED
PATENTED
CERTIFIED
SAFE



Notified Body n° 0426

est report n° PPE125RP2602 Page 4 of 4

The values of breathing resistance (inhalation and exhalation) were recorded continuously and are shown below

Results

Id.sample	Inhalation resistance ΔP (mbar)		Exhalation resistance ΔP (mbar)	
	25 x 2 Vmin	25 x 2,5 l/min	25 x 2 l/min	25 x 2,5 l/min
#A1	6,4	9,2	5,8	7,2
#A2	6,3	8,3	6,0	7,4
#A3	5,7	8,1	5,4	7,2







ARIA 19 FIGHTS AGAINST COVID 19

The spread of COVID-19 pandemic sadly led to PPE shortages all around the world, causing a tremendous challenge to the U.S. healthcare system.

Since 1980, OceanReefGroup has been a worldwide leader in full face masks and the diversification of our expertise allowed us to design and produce a mask with an innovative breathing filtration system and a full coverage face shield.

THE PROJECT

In the beginning of March 2020, we deposited an international Patent through our company Mestel Safety, under the brand OCEANREEF manufactures full face masks purposes and through the brand Safety Protection Equipment, manufactures gas masks.



DEVELOPMENT OF THE SAFETY STANDARD BEHIND OUR PRODUCT



Upper volume (RED) and Lower volume – technically called Orinasal Pocket (BLUE)

The FFM shares this same design, which means that the breathing cycle of both products perform in the same manner.

BREATHING - CO₂



- 1. Inhaled air goes from outside to inside the mask from the filter
- 2. Inhaled air is transferred from the upper volume (where it fogs visibility through the mask) to the lower volume, via the one-way valves in the orinasal pocket, reaching the inside of the orinasal pocket.
- 3. The air is, then, inhaled by the user.
- 4. Exhaled air cannot move back into the upper volume, because the orinasal pocket seals the nose and mouth around the cheeks. The valves are one-way valves and do not open upwards.
- 5. The exhaled air moves in the only direction it can: through an exhalation one-way valve.

This breathing cycle is called one-way breathing circulation. The benefit is to separate fresh air, containing only 0.4% CO2, from used, exhaled air with 4% CO2 concentration.

FEATURES



- Easy use APA and RD40 adaptors are easy to assembly on the Aria 19 mask, also when it is weared.
- The strap can be quickly adjusted and fits well almost any face shape, allowing a comfortable wearing and removing.
- The voice can be heard clearly throughout the Aria 19 mask thus allowing a good communication.
- The high-transparency visor allows a perfect vision (wide optical field) with a slight optical distortion in few areas.



VARIOUS COMPATIBLE FILTERS



CO₂ TESTED



REUSABLE



FOG FREE

CORPORATE SOCIAL RESPONSABILITY

Thanks to Aria19 KIT we were able to give a small contribution to help thousands of health care in the right against coronavirus.



ORGANISATION CHART

MEDICAL DIVISION



Marketing Manager

Luca Gamberini

I.gamberini@oceanreefgroup.com



Business Development Manager

Piergiuseppe Donato

p.donato@oceanreefgroup.com

Customer Care

Europe, Middle East, Asia, Oceania



• Samantha Moresco

s.moresco@oceanreefgroup.com



Federico Giunto

f.giunto@oceanreefgroup.com



Sales & Information

Europe, Middle East, Asia, Oceania Margherita Lunetta

megl@oceanreefgroup.com



North & South America
Sandra Campos-Montalvo
operator@oceanreefgroup.com



Customer Service/ Tech Support

Europe, Middle East, Asia, Oceania Fabio Porcile

f.porcile@oceanreefgroup.com



North & South America
Jesse Negrete
customerserviceusa@oceanreefgroup.com



North & South America Eleonora Fontana e.fontana@oceanreefgroup.com