

#### Water Distiller

Use mineral-freewate r to protect your sterilizer and instruments at a fraction of the cost of bottled distilled water. It passes all safety tests, complying with all electromagnetic compatibility standard (2014/30/EU) and low voltage standard (2014/35/EU) requirements.



#### Steam Sterilizer

Sterilizers and accessories that help dentists and their assistants to always ensure the highest levels of safety and hygiene in their practice



Manufacturer & After service Ningbo JOIDENT Electronics Technology Co., Ltd. No. 111 JinGuNan Road, Yinzhou, Ningbo 315100, P.R.China

t +86-574-8732 8073 www.joident.com f +86-574-8732 9073 info@joident.com




The sealing machine has no complicated operation, simple Operation, reasonable design, safe and reliable size. Moderate quality, light weight and high efficiency




### **Sealing Machine**

80+ Countries User 12 months warranty 5+ Patent technology 1800+ days Aging test

Recorded in December 2020.







# MODELING

Combining the experience of modern art and minimalist shape design, it is seamless from inside to outside. Our design size can be used for any ocation, so that you can experience the charm of modern art.
















**OPERATION** 

# **LEADING EDGE**



Bi-directional

Pros and cons are available, ergonomic design



Safety slip

The sealing machine composed of a front indicator and a sound reminer




**Automatic Alarm** 

Hidden sliding blade to protect the user from injury



# **ANTI-ROLL AND ANTI-SCORCH**

The sealing quality (12mm-wide sealing seams) is intact after 3 seconds of pressure test with no risk of burning



#### **ONE PREHEATING CAN BE USED CONTINUOUSLY**









The sealing temperature and pressure can be adjusted automatically by the intelligent sealing system that the machine can still keep high performance sealing under continuous use




Flexible adjustment according

To the width of sterilization bag



large selection of Roll width at the same time (3 rolls), with high efficiency




The maximum bearing width is up to 300mm



ROLL HOLDER ALL TYPES OF ROLL



### **Exquisite Accessories**



Microcomputer Control system



Condensing system

Condensation is provided


The optional function of the system


Optional




Instrument tray
Optional

#### **Product Size**







#### **Technical Data**

Voltage frequency: 220V ,110V/50HZ,60HZ

Consumption Power: 100W Preheating time 2-3minutes

Packing width: 12mm longest seal: 30cm Overheat protection: Yes

Packaging size: 350\*350\*535CM

Gross weight: 7.5KGS

Phase number

Water proof grade I/IPXO

Fuse specification: 2-5x 20 T8A

· Indoor use

· Altitude less than 3000m Service Environment

· Ambient temperature: 5 to 40 degrees

· Relative humidity: 85%RH · Maximum fluctuation:+-10%

### **Conformity**

Low Voltage Directive	2014/35/EU
Electromagnetic Compatibility Directive	2014/30/EU
Machinery Safety Directive	2006/42/EC

## **MEDICAL STERILIZATION BAG**

Sterilisation is one of the most important measures involved in rigorous dental practice hygiene. Hygofol sterilised packaging offers maximum reliability with its multi-layer laminated foil and highly tear-resistant sterilisation paper. In combination with the JOIDENT sealing machine process sealing device, this product can be used to create easily peelable 12 mm-wide sealing seams. This enables you to meet all requirements in accordance with EN ISO 11607 and EN 868.



200m/roll	Weight	Order number
5cm	1.23kg	00502
7.5cm	1.85kg	07752
10cm	2.5kg	01002
15cm	3.7kg	01502
20cm	5.3kg	02002
25cm	6.25kg	02502
30cm	7.4kg	03002