



ticare TISSUE CARE
PHILOSOPHY

THE RESULT OF A UNIQUE
MANUFACTURING PROCESS
AND QUALITY CONTROL:

Manufacturing
tolerances up to
7 microns for a
perfect
adjustment.

11 degrees conical
connection for a
bone respectful
load distribution
and perfect
seal.



gapZero
ticare

THE IMPLANT FREE OF
BACTERIAL MICROLEAKAGE

The **RESULT:**
the implant
free of bacterial
microleakage that
reduces the risk of
periimplantitis and
maintains the marginal
bone level.

100% verification
to guarantee the
highest implant and
abutment quality
standards.





ticare

TISSUE CARE
PHILOSOPHY

It is two years now since the publication of the Scientific Study that proves ticare is

THE IMPLANT FREE OF BACTERIAL MICROLEAKAGE

BACTERIAL MICROLEAKAGE AT THE ABUTMENT-IMPLANT INTERFACE, IN VITRO STUDY. Larrucea C. et al. Clinical Implant Dentistry and Related Research 2018; 1-8.



Scan the code to see the explanatory video of the study
ticareimplants.com/gapZero



ticareimplants.com

@ticareimplants

+34 983 309 602

info@ticareimplants.com



Groups with 20 N and 30 N torque applied showed no bacterial leakage at any time during the experiment.

