



Kagayaki® Enforce Pin

3-step polishing system made of Urethane for contouring, smoothing and finishing of composite, compomer and glass-ionomer restorations

40 pcs/box. Assortment 30 pcs/box.

3000-5000 RPM	for fast removal and contouring			for smoothing and silky gloss			for achieving an ultra high shine		
grit size	① coarse (125μ)			② mid-coarse (70 μ)			③ fine (32 μ)		
colour	white			pink			yellow		
REF	EP 125-1	EP 125-2	EP 125-3	EP 70-1	EP 70-2	EP 70-3	EP 32-1	EP 32-2	EP 32-3
	assortment EP 125			assortment EP 70			assortment EP 32		

Kagayaki® Enforce

3-step polishing system made of Urethane for contouring, smoothing and finishing of composite, compomer and glass-ionomer restorations. Delivered with mandrel.

25/50 pcs/bag.

3000-5000 RPM	for fast removal and contouring			for smoothing and silky gloss			for achieving an ultra high shine		
grit size	① coarse (125μ)			② mid-coarse (70 μ)			③ fine (32 μ)		
colour	white			pink			yellow		
REF	E 125-1	E 125-2	E 125-3	E 70-1	E 70-2	E 70-3	E 32-1	E 32-2	E 32-3
	assortment E 125			assortment E 70			assortment E 32		



Kagayaki® Ensmart Pin

3-step polishing system made of silicone for contouring, smoothing and finishing of composite restorations.

30 pcs/box.



3000-5000 RPM	for fast removal and contouring			for smoothing and silky gloss			for achieving an ultra high shine		
	① coarse (125μ)			② mid-coarse (70 μ)			③ fine (32 μ)		
	white			green			pink		
	REF	ENP 125-1	ENP 125-2	ENP 125-3	ENP 70-1	ENP 70-2	ENP 70-3	ENP 32-1	ENP 32-2
	assortment ENP 125			assortment ENP 70			assortment ENP 32		

Kagayaki® Ensmart

3-step polishing system made of silicone for contouring, smoothing and finishing of composite restorations. Delivered with mandrel.

25/50 pcs/bag.



3000-5000 RPM	for fast removal and contouring			for smoothing and silky gloss			for achieving an ultra high shine		
	① coarse (125μ)			② mid-coarse (70 μ)			③ fine (32 μ)		
	white			green			pink		
	REF	EN 125-1	EN 125-2	EN 125-3	EN 70-1	EN 70-2	EN 70-3	EN 32-1	EN 32-2
	assortment EN 125			assortment EN 70			assortment EN 32		



Kagayaki® Ensmart Pin

3-step polishing system made of silicone for contouring, smoothing and finishing of composite restorations. Stainless steel pin.

30 pcs/box.



3000-5000 RPM	for fast removal and contouring			for smoothing and silky gloss			for achieving an ultra high shine		
	① coarse (125μ)			② mid-coarse (70 μ)			③ fine (32 μ)		
grit size	white			green			pink		
colour	ENPS 125-1	ENPS 125-2	ENPS 125-3	ENPS 70-1	ENPS 70-2	ENPS 70-3	ENPS 32-1	ENPS 32-2	ENPS 32-3
REF									



Kagayaki® Export Single step polishers for composite

More pressure for initial finish, less pressure for smoothing.

30 pcs/box.



3000-5000 RPM	polishing and finishing		
grit size	coarse (125μ)		
colour	whitish		
REF	1003-E	1002-E	1001-E



Kagayaki® Diamond Sun

Diamond impregnated polisher for high surface luster of composite restorations.

30 pcs/ box.




<p>3000-5000 RPM</p> <hr/> <p>grit size</p> <hr/> <p>colour</p> <hr/> <p>REF</p>	for achieving a smooth surface			for achieving an ultra high shine, enhancing esthetic properties of the filling		
	① fine (30 μ)			② extra-fine (10 μ)		
	beige			yellow		
	DS 30-1	DS 30-2	DS 30-3	DS 10-1	DS 10-2	DS 10-3
assortment DS 30			assortment DS 10			



Kagayaki® RoundFlex

Polishing and finishing complete set. Delivered with 1 short and 1 long mandrel

50 pcs./box 100 pcs./box



<p>20000 RPM</p> <hr/> <p>grit size</p> <hr/> <p>colour</p> <hr/> <p>dia.</p> <hr/> <p>REF</p>	contouring		smoothing		polishing		high gloss	
	coarse		medium		fine		super fine	
	brown		red		orange		yellow	
	Ø12,6 MM	Ø9,6 MM	Ø12,6 MM	Ø9,6 MM	Ø12,6 MM	Ø9,6 MM	Ø12,6 MM	Ø9,6 MM
2028C	2027C	2026M	2025M	2024F	2023F	2022SF	2021SF	
Assortment Ø 12.6 mm: 100 pcs(25 pcs each type+2 free mandrels)						2030/40		
Assortment Ø 9.6 mm: 100 pcs(25 pcs each type+2 free mandrels)						2029/40		



Kagayaki® Mirā

Polishing paste:
Fine polishing of fillings.
Removal of food and tobacco stains.

- mint (REF 3016) 45g

45 g. jar



Kagayaki® Mirā Zr (with Zirconium)

Polishing paste:
Fine polishing of fillings.
Removal of dental tartar deposits, food and tobacco stains.

- mint (REF 3012) 50g
- mint (REF 3013) 120g

50 g. jar
120 g./tube



Kagayaki® Henka

One step prophy-paste offering both an efficient cleaning and polishing effect.
Due to special abrasive particles Henka provides a cleaning action during the first seconds of application with a gentle yet powerful polishing action.

Available in several flavours:

- apple (REF 3021)
- mint (REF 3022)

115 g./tube



Kagayaki® Kurin One. Paste for cleaning

Paste for quick, easy and effective natural teeth cleaning. Ideal for metallic and composites filling cleaning as well.

- anise (REF 3031)


50 g./ tube

Kagayaki® Kurin Two. Paste for polishing

The ideal product for natural teeth and metallic fillings super high-gloss polishing.

- anise (REF 3032)

50 g./ tube




Cups for restoration polishing in conjunction with Kagayaki® polishing paste.

Polishing cups to use in conjunction with Kagayaki paste for restoration polishing. The polishing paste is absorbed into the pores of a sponge and stays there during the polishing procedure which enables fast polishing and decreases the consumption of polishing paste. Delivered with mandrel.

REF 1020/4S.

25 pcs/ bag.

KAGAYAKI® Endo




Ø 0.8	108	208	308
REF	108W	208Y	308R

Kagayaki® Burnout posts YuniPin

Burn-out plastic posts.

100 pcs/bag.



Kagayaki® Silicone stops. Autoclavable

Silicone-Stop Dispenser.

colour	yellow	blue	red
REF	SDY	SDB	SDR

100 pcs/bag.