

Carefix[®]

made by



A Superior Barrier Mask for civilian use offering over 90% protection.

We are proud to introduce Tytex's Barrier Masks.

Carefix[®]

Tytex A/S, a 50 year-old Danish company, is a global leader in medical textiles and devices.

www.tytex.com

Tytex products meet the highest quality standards which are approved by health care authorities around the world.

Benefits

Medical grade fabrics and textiles contribute to the quality of the product and to its superior level of protection when compared to homemade masks and commercially made paper and cloth masks.

- **Filtration Efficiency:** Protects wearer against airborne viruses and other contaminants.
Protects wearer against hand-to-face contamination.
Protects others if the wearer is asymptomatic.
- **Test Results:** 97% Effective Pre-Wash; 93% after 25 Washes.
- **Reusable:** Washable up to 25 times.
- **Excellent Value:** Compared to other single use disposable masks.
- **Comfort:** Comfortable to use over prolonged time periods.
Easy to put on and remove; no elastic bands around ears.
- **Convenience:** May be worn around neck when not in use.

Features

Carefix[®] Barrier masks are used primarily for non-medical everyday civilian use. They are not a substitute for N95 or surgical masks but may be used on top of N95/Surgical masks to extend their useful life.

- Available Sizes: Small / Medium and Medium/Large.
- Colors: White, Black and Tan.
- Packaging: Retail: 5 Masks per Bag (in cases of 400 masks).
Bulk: 100 Masks per Bag (in cases of 400 masks).
- Composition: Medical Grade: Polyamide 94% & Elastane 6%.
- Washing Instructions: Washable at normal washing temperatures.



Carefix®

Carefix[®]





Carefix[®]



Customer logos can be woven into the mask at no additional cost.



Inside View (Label Reversed)



Inside View (Label Reversed)

Instructions for use

Make sure to choose the right size. The mask must sit tight around your mouth and nose.

Carefix®

Application

- Wash your hands thoroughly prior to application of the mask.
- Apply the mask from the top of the head with the wide part toward the nose and the narrow part toward the neck.
- With both hands, pull the mask down around the head, making sure to cover the nose and mouth.
- Arrange the mask just under the ears and around the neck. Be careful not to touch mouth, nose or eyes when applying the mask.



Removal

- Wash your hands thoroughly before the mask is removed.
- Remove the mask by lifting the narrow part covering the neck over the head and then remove the mask from mouth and nose, while still holding on to the narrow part. Be careful not to touch mouth, nose or eyes when removing the mask.



TEST RESULTS

Standard: EN13274-7

Test Item: Filtration of Aerosols

Sample: A

Result

Requirement

Before Wash:

97.0%

UNS 1* \geq 90.0%

After 25 Washes:

93.0%

* UNS 1: For people who have frequent and regular contact with the public or other people, such as policemen or cashiers.

Contact information

For additional information regarding:

- Minimum Order Quantity.
- Delivery Lead Times.
- Product Availability.
- Pricing.
- Any questions you might have.

Contact: JD Honigberg International, Inc.

Northbrook IL 60062 USA

Email: medical@jdhintl.com

Phone: +1-847-412-0200

www.jdhmedical.com