

Zirconia One-Stop Solution



Aconia TT Multilayer / TT White



- Extreme incisal translucency
- Singles to three-unit anterior bridges
- 16 VITA Classical Shades

TRANSLUCENCY
46-49%

STRENGTH
600-900Mpa

Aconia SHT Multilayer / SHT Pre-shaded / SHT White



- Translucent with reliable durability
- Wide indications for full contour
- 16 VITA Classical Shades and 3 Bleach Shades

TRANSLUCENCY
43-46%

STRENGTH
900-1100Mpa

Aconia ST Pre-shaded / ST White



- Strong with good millability
- Wide indication from coping to long- span bridge
- 16 VITA Classical Shades

TRANSLUCENCY
43%

STRENGTH
1250Mpa

Aconia HT+ White



- Best option for framework & abutment
- Economical solution

TRANSLUCENCY
42%

STRENGTH
1400Mpa

Coloring Liquid & Other Consumables




Aconia Coloring Liquid

Standard



- Precisely matched to VITA 16 color system
- Well-suited for both dipping and brushing (paint-on) methods



Master

- Precisely matched to VITA 16 color system
- Ideal results applied on Aconia white blanks by Aconia Coloring Technology

Artist



- Reproduce rare and special colors with esthetic art effect, natural and lifelike
- Start coloring process without preparing and mixing

Other Consumables



PMMA

Easy to polish
Various color available
Long-term shade stability



WAX

Easy to mill
Burn out completely without residue
Compatible for all male-molds oral tissue



MILLING BURS

Excellent durability
Diamond & Uncoated
Main systems available



Chengdu Besmile Biotechnology Co., Ltd.

Tel: +86-28-85317108 | Fax: +86-28-85217510 | www.bsmdental.com

No.1666, 2nd section of Chenglong Road, Chengdu Economic and Technological Development Zone, Sichuan, China

